

Data collection and data reduction

Which projects are involved?

- *DNA* - automation of whole process from crystal to $|F_{\text{obs}}|$ and beyond...
- *BEST* - analytical strategy program
- *Mosflm* - improvements to data integration
- *Mosflm* - *ccp4i* task
- *Scala* - improvements to scaling

DNA - what's it supposed to do?

- Automate data collection and processing at synchrotron beamlines
- Use & develop existing software to automate data processing
- Write new software to coordinate DC and DP

<http://www.dna.ac.uk>



DNA - who's involved (where)?

BEST & Mosflm


- *BEST* is a strategy program which uses more information than the current strategy option, and provides more detail for the optimum data collection
- currently being;
 - modified to use output from *Mosflm*
 - incorporated into *DNA* using output from *Mosflm*

Sasha Popov and Gleb Bourenkov (*EMBL, Hamburg*)

Mosflm - who does what?

- Andrew Leslie - coordinates whole thing
- Harry Powell (Andrew Leslie) - bug fixing, maintenance, development, user support, *etc, etc.*
- Geoff Battye (Harry Powell) - new GUI

complete downloads of Mosflm v6.2.2 & v6.2.3 by month



Mosflm - recent past & short term future

- major new developments
 - new GUI
 - spot positions in postrefinement
- minor new features
 - ice-ring exclusion in autoindexing
 - performance improvement (buffered *i/o* for SPOTOD)
 - *o/p* modified for *BEST* & *DNA*
- improvements
 - mosaicity estimation
 - autoindexing
- maintenance
 - bug fixing, additional detectors (automatic recognition)...
 - user support, better diagnostics

New *Mosflm* GUI

- Designing a new GUI for *Mosflm* with the aim of providing an attractive and straightforward interface to data integration in a step-by-step approach.
- communicates with *Mosflm* via *TCP/IP* sockets using a custom XML

Geoff Battye (*MRC-LMB, Cambridge*)

New GUI - timescale

- alpha test version (locally in LMB) - end June
- beta release (developers) - end September
- beta release (public) - before Christmas

☒ Edits allowed

Processing params

a : 0.00
 b : 0.00
 c : 0.00
 alpha : 0.00
 beta : 0.00
 gamma : 0.00
 PsiX : 0.00
 PsiY : 0.00
 PsiZ : 0.00
 Mosaic : 0.300
 Divh : 0.000
 Divv : 0.000
 Lambda : 1.542
 Distance: 194.20
 Beam X : 149.78
 Y : 151.45
 CCOMEGA : 0.000
 ROFF : 0.00
 TOFF : 0.00
 YSCAL : 1.0000
 Pick area: X: 11
 Y: 11
 Int threshold: 20
 Vector scale 1
 Two theta 0.00
 Resolution 0.00
 SPOT SEARCH
 Threshold 5.00
 Rmin 15.00
 Rmax 135.00
 X offset 0.00
 Y offset 0.00
 Min X size 0.50
 Max X size 2.00
 Min Y size 0.50
 Max Y size 2.00
 Min no of pix 6
 X splitting 0.30
 Y splitting 0.30
 AUTOINDEXING
 Min I/sig(I): 20

Prompts ☒ OnUpdate display: ☒ NoAfter refinement ☒ NoAfter integration ☒ NoTimeout mode ☒ Off☒ Select item

Main menu

Read image
 Find spots
 Edit spots
 Clear spots
 Select images
 Autoindex
 Estimate mosaicity
 Predict
 Clear prediction
 Adjust
 Refine cell
 Integrate
 Strategy
 Keyword input
 Find hkl
 Pick
 Measure cell
Circles
 Beam / backstop
 Save/Exit


☐


Output

Pixel X,Y 2001 0
 XC,YC mm 300.0 0.0
 Resolution 0.00
 Indices 0 0 0
 F Phi 0.00 width 0.00
 Intensity 0
 Sigma 0
 Spacing A 0.000
 Average 0.0
 Rms 0.0
 Number 0
 Zoomfactor 0
 Circle resolution A
 0.0 0.0 0.0 0.0
 Phi 0.00 0.20
 Missets ThetaX,Y,Z
 0.00 0.00 0.00

☒ lys_fine_002.pck

Min 1 Max 253 Cursor position

Overlay ☒ on Contrast Colour ☒ Black on whi Mag ☒ x4 



Mosflm ccp4i task

- designed only for integrating images - no functionality is intended for indexing or postrefinement as these are best performed under supervision

Pete Briggs (*CCP4, Daresbury*)

Run MOSFLM in batch mode

Help

*Run MOSFLM in batch mode to integrate diffraction images*Job title Load parameters from command file Browse View

Load parameters from command file

Specify where the images are located☒ Use current working directory /gr0/harry/lys_fineImage file name template: Matrix file trypsin Browse ViewMTZ out trypsin Browse View*Dataset and Harvesting Parameters*☐ Create harvest file in project harvesting directoryCrystal belonging to Project trypsinDataset name *Images to Integrate*Edit list Add processing run*Crystal Parameters*Space group Cell a b c alpha beta gamma Mosaicity Wavelength of radiation used in the experiment Angstrom*Detector Parameters*☐ Set detector type to Mar Research IP and over-ride default determined by MOSFLMRun Save or Restore

Close

Scala - progress


- rewrite using OOP methods
- unmerged reflection class which can be used for scoring all possible point groups (initially in a standalone application, later to be incorporated in Scala)

Phil Evans (*MRC-LMB, Cambridge*)


Scala_MPI - parallel version

- Improved performance by parallelizing the rate determining steps.
- Main bottleneck is the construction of a parameter matrix from the reflection data which is used to scale and merge the data; it is constructed by looping over all the observations in the MTZ file.
- Each available processor takes its own section of the MTZ file and compiles its own parameter matrix. These are then collected and summed into a single matrix.
- Code parallelised using MPI for use on distributed memory clusters.

Ronan Keegan (*eHTPX*, *Daresbury*)



Scala_mpi: speedup of parallel code on up to 4 processors



that's all folks!