



# VMXm: Micro/nanofocus protein crystallography beamline at Diamond


Anna Warren

[anna.warren@diamond.ac.uk](mailto:anna.warren@diamond.ac.uk)



# Introduction

- Target proteins are getting more complex, leading to smaller and more disordered crystals
- Can design ideal experiment to give optimal data quality



# Current Limits

- 2.2 Å data can be collected from 1  $\mu\text{m}^3$  crystals (~700 well diffracting crystals)
- 3 Å data can be collected from 5  $\mu\text{m}^3$  membrane protein crystals grown in LCP (~35 crystals, grid scanned first for centring)
- From theoretical calculations a complete 2 Å dataset can be collected from a single 1  $\mu\text{m}^3$  lysozyme crystal (Holton and Frankel, 2010)
- Discrepancies between theory and experiment

# Current Limits

- Dose tolerance of samples cannot be changed – Henderson/Garman limit fixed
- Reduce dose on sample to measure given data quality:
  - Reduce experimental background
  - Cleaner sample mounting
  - Improve analysis for weak and multicrystal data
  - Record rotation data to improve data quality
  - Visualization of micron and sub-micron crystals
- Take advantage of photoelectron escape



# VMXm Aims

- Improve signal to noise by reducing background:
  - Sample environment under vacuum
  - Crystals mounted with minimal high quality rotation diffraction data
  - Reduce beamsizes to match that of the crystal where both **size** and **number** of crystals are limited
- Standard rotation data collection on samples down to 500 nm
  - Alignment without the need for X-ray raster scanning
- Optimise sample alignment, sample cooling and data analysis for micron and sub-micron crystals
- Data collections using minimal amounts of sample

# VMXm Specifications

- 6 – 28 keV energy range
- 0.3 – 10  $\mu\text{m}$  (v) & 0.5 – 5  $\mu\text{m}$  (h)
- Flux
  - $\sim 10^{12}$  ph/s in  $0.3 \times 5 \mu\text{m}$  (v x h)
  - $> 10^{11}$  ph/s  $0.3 \times 0.5 \mu\text{m}$  (v x h)
- Interchangeable detectors:
  - Pilatus3 6M Si
  - Eiger2 X 9M CdTe

# Sample Environment

- *In vacuo* sample environment - reduce X-ray background to a minimum
- Cryo-stage - preserve sample *in vacuo*
- Standard on-axis optical microscope + SEM for sample visualization and alignment
- High stability goniometry - permit rotation data to be measured from micro and nanocrystals
- Serial data collection also possible



# Sample Preparation






# Sample Preparation


Ideal sample mounting:

- Reduce background to increase signal-noise:
  - Air path around sample
  - Sample mount
  - Solvent surrounding crystal
- Multiple crystals per mount
- Cryo-EM style sample preparation
  - Grids
  - Blotting
  - Plunge cooling
- SEM for sample characterisation




# Cryo-EM Grids



3 mm






Support film





Grid











# Sample Blotting



# Sample Characterisation





# Sample Mounting



# Sample Preparation Lab








# Beamline





# Current Status





Cryo gripper

Fluorescence detector

Sample holder

OAV port

Beamstop XY

Sample XZ stages

Goniometer

Isolated cold-stage and gripper

SEM

SED






# Beamline Results




# Background Measurements

$< 0.1$  counts per pixel per second



Atmosphere





Vacuum



# I24




# VMXm





# CPV Ld14 – Cytoplasmic polyhedrosis virus

- Spacegroup I23
- Unit cell  $a=b=c=103 \text{ \AA}$
- Approximate crystal size  $2 - 4 \text{ \mu m}$
- Data collected at 21.3keV – exploit photoelectron escape
- Eiger2 X CdTe 9M





# CPV Ld14 – Cytoplasmic polyhedrosis virus






# CPV Ld14 – Cytoplasmic polyhedrosis virus

	Ji et al., I24	VMXm
Detector	Pilatus 6M Si	Eiger2 X CdTe 9M
Number of crystals	20	1
Energy keV	12.458 keV	21.3 keV
Resolution	72.5–1.91 (2.02-1.91)	72.9-1.30 (1.32-1.30)
Unique Reflection	12952(1045)	40616(1350)
Completeness (%)	92.3(52.1)	90.7(62.2)
R <sub>merge</sub>	0.199(0.327)	0.284(0.881)
I/ $\sigma$ I	8.4(1.8)	3.8(0.5)
CC1/2	??	0.933(0.407)
Beamsize	6 x 6 $\mu$ m slit	3.6 x 3.6 $\mu$ m



# CPV Ld14 – Cytoplasmic polyhedrosis virus

Stat	Value
Resolution	72.94-1.38
N. Reflections all/free	37547/1920
R/Rfree	0.136/0.193
RMS dev	
Bonds	0.0141
Angles	1.806



# Membrane protein/LCP samples

## Challenges of LCP:

- Viscous, hydrophobic - difficult to handle.
- Support films on cryoTEM grids are fragile.
- Generates significant background signal in X-ray diffraction experiments.


## Changing the phase: Cubic → Sponge:

- Sponge has reduced viscosity - it is 'pipettable'
- Maintains a hydrophobic environment.
- Smaller molecular weight precipitants can be used as 'spongifying' agents:
  - MPD, PEG200, PEG400



# FFAR1


- GPCR membrane protein
- Grown in classic glass plate sandwich - generates crystals in variety of sizes
- Crystallisation solution including 10% MPD used to wash the bolus out from under the glass.
- Applied directly to grid mounted in plunge freezer as per soluble sample.



# FFAR1

- Data collected from 7 crystals at 21.3 or 19.5 KeV and merged.
- Completeness potentially affected by orientation and C2 spacegroup.

Number of crystals	7
Energy (keV)	21.3 or 19.5
Resolution (Å)	100.98 – 2.66 (2.71-2.66)
Spacegroup	C2
Unique Reflections	14118(737)
Completeness (%)	76.1(78.6)
R <sub>merge</sub>	0.558(2.842)
I/σI	3.3(0.5)
CC1/2	0.913(0.155)



Resolution (Å)	100.98-2.66
No. reflections all/free	14152/680
R-factor/R-free	0.248/0.293
RMS Deviations	
Bonds (Å)	0.0054
Angles (°)	1.460



# Final Remarks

- VMXm aims to record high quality diffraction data where both size and number of crystals are limited
- Exploiting low background, sample mounting, high energy and photoelectron escape to optimize experiment
- Successfully demonstrated the VMXm principal of operation
  - Rotation data recorded from crystals a few microns in size
  - Crystals mounted on cryoEM grids in vacuum at cryo-temperatures

# Acknowledgements

Graham Duller  
Andy Stallwood  
Richard Parr  
Tony Gardner  
Dan Magrath  
Marcell Nagy

Leo Cahill  
Richard Littlewood  
Kevin Wilkinson  
Alan Nash  
Steve Buxcey  
Paul Benson's IFM Team

Andy Foster  
James O'Hea  
Brian Nutter  
Ben Kemp

Victoria Lawson  
Neil Smith  
Paul Coe

Mark Lunnon  
Dave Butler  
Adam Taylor  
Adam Prescott  
Leon Adams  
Thomas Hartrampf

Simon Lay  
Hugo Shiers  
Guenther Rehm  
Chris Bloomer  
Sandira Gayadeen

David Laundry  
Lucia Alianelli  
John Sutter  
Simon Alcock  
Kawal Sahwney

Sonia Moon

## **VMXm**

**Gwyndaf Evans**

Jose Trincao  
Adam Crawshaw  
Emma Beale

## **STRUBI, Oxford**

David Stuart  
Geoff Sutton

## **Membrane Protein Lab**

Andrew Quigley  
Harish Cheruvara  
James Birch



# We need you!

- Do you have suitable crystals with dimensions  $<10\text{ }\mu\text{m}$ ?
- Beamline staff will assist with sample preparation and data collection
- Please contact us to discuss:
  - [anna.warren@diamond.ac.uk](mailto:anna.warren@diamond.ac.uk)