


Data Processing with XDS

Kay Diederichs



Outline

1) Introduction to XDS

2) Data processing is the crucial link between your experiment and your structure

- * garbage in – garbage out!
- * for important data, do not rely on automatic data processing alone
- * manual checks are easily performed with XDSGUI
- * try to optimize data processing – this converts noise to signal, and may enable structure solution, and/or improve refinement
- * understand statistics (next talk)

3) Process your data (afternoon, and during the week)

Tutors: Tiankun Zhou, Ralf Flaig, Clemens Vornrhein, me

The *XDS* program suite

- Original author: Wolfgang Kabsch (Max-Planck-Institute Heidelberg)
- Since ~1986
- I joined 2007



The XDS+ programs

- **XDS**: the main program - indexing, integrating, scaling, statistics
- **XSCALE**: scale several XDS intensity data sets together; zero-dose extrapolation; statistics
- XDS CONV: convert to MTZ / SHELX /... format (AIMLESS and CTRUNCATE are not needed!)

Programs independent from the XDS distribution:

- XDS-Viewer: inspect diagnostic images written by XDS, or (single) data frames (open source). *adxv* or *dials.image_viewer* may be used instead.
- XDSSTAT: additional statistics
- **XDSGUI**: graphical user interface for XDS, SHELX C/D/E, ARCIMBOLDO (open source)
- **XDSCC12**: (XDS) which frames are bad?
(XSCALE) which data sets to re-index and merge?

Sources of information

- XDS main website: <https://xds.mr.mpg.de> ; download for Linux (WSL), Mac. complete, accurate, up-to-date documentation
- XDSwiki: <http://strucbio.biologie.uni-konstanz.de/xdswiki>
Installation; data sets; documentation; download; links to e.g. Matthew J. Whitley's excellent tutorial given at CSHL 2018
- CCP4 bulletin board
- SBGrid talk (May 2020) at <https://www.youtube.com/watch?v=3WU9NrILECo>
- Making a difference in multi-data-set crystallography : simple and deterministic data-scaling/selection methods. Assmann, G.M., Wang, M., Diederichs, K. (2020) Acta Cryst D76, 636 (serial crystallography, XDSCC12)

Automatic processing with XDS

- beamline software (provides **XDS.INP**)
- scripts: **xia2** (CCP4), **autoPROC** (Globalphasing), **generate_XDS.INP** (XDSwiki), **fast_dp** (Diamond), *xdsme* (Soleil), *autoxds* (SSRL), *autoprocess* (CMCF), ...
- CCP4: *pointless*, *xdsconv* (type CCP4_I+F, or CCP4, or CCP4_I, or CCP4_F)
- SHELX: *shelxc* reads XDS_ASCII.HKL

Principle of XDS processing

- There is one JOB= line in **XDS.INP** which specifies a list of tasks:

JOB= XYCORR INIT COLSPOT IDXREF DEFPIX INTEGRATE CORRECT

- data reduction is divided into tasks in a **modular** way
- information storage/exchange/flow between tasks by data files which may be inspected/analyzed
- each task needs the result from the previous tasks
- fine-tuning of a task does *not* require previous tasks to be repeated


– each task writes its output file **<TASK>.LP**

The tasks are ...

- XYCORR : write positional correction files
(**X-CORRECTIONS.cbf**, **Y-CORRECTIONS.cbf**)
- INIT : find background pixels (defaults usually OK)
- COLSPOT: find reflection positions
- IDXREF : "index" reflections; user may supply/choose spacegroup
- XPLAN [not required] : strategy for data collection
- DEFPIX : mask shadows on detector (use *XDSGUI!*)
- INTEGRATE : evaluates intensities on all frames, writes **INTEGRATE.HKL** and **FRAME.cbf**
- CORRECT : **scales**, rejects outliers, statistics, writes scaled, unmerged **XDS_ASCII.HKL** (and other files)

Information flow

NAME_TEMPLA OSCILLATION_RORGX DATA_RANGE
TE_OF_DATA_F ANGE ORGY
RAMES SEPMIN DETECTOR_DISTANCE
DETECTOR STRONG_PIXEL X_RAY_WAVELENGTH
SPACE_GROUP_NUMBER



Example **XDS.INP**

```
JOB= XYCORR INIT COLSPOT IDXREF DEFPIX INTEGRATE CORRECT
ORGX=1546 ORGY=1552      !Detector origin (pixels); e.g. NX/2 NY/2
DETECTOR_DISTANCE=180    ! (mm)
OSCILLATION_RANGE=0.50   !degrees (>0)
X-RAY_WAVELENGTH=0.980243 !Angstroem
NAME_TEMPLATE_OF_DATA_FRAMES=frms/wga2-27_1_???.img
DATA_RANGE=1 360         !Numbers of first and last data image collected
BACKGROUND_RANGE=1 10    !Numbers of first and last data image for background
SPACE_GROUP_NUMBER= 19   !0 for unknown crystals; cell constants are ignored.
UNIT_CELL_CONSTANTS= 44.4 86.4 104.5 90 90 90 ! not required if spgr=0
REFINE(IDXREF)=BEAM AXIS ORIENTATION CELL DISTANCE
REFINE(INTEGRATE)=DISTANCE BEAM ORIENTATION CELL ! AXIS
ROTATION_AXIS= 1.0 0.0 0.0
INCIDENT_BEAM_DIRECTION=0.0 0.0 1.0
FRACTION_OF_POLARIZATION=0.99                                ! SLS X06SA
POLARIZATION_PLANE_NORMAL= 0.0 1.0 0.0
DETECTOR=CCDCHESS      MINIMUM_VALID_PIXEL_VALUE=1          OVERLOAD=65000
DIRECTION_OF_DETECTOR_X-AXIS= 1.0 0.0 0.0
DIRECTION_OF_DETECTOR_Y-AXIS= 0.0 1.0 0.0
VALUE_RANGE_FOR_TRUSTED_DETECTOR_PIXELS= 7000 30000 !Used by DEFPIX
                                           !for excluding shaded parts of the detector.
INCLUDE_RESOLUTION_RANGE=50.0 1.3 !Angstroem; used by DEFPIX,INTEGRATE,CORRECT
```

Bold keyword/parameter pairs are required. Complete documentation at
xds.mr.mpg.de/html_doc/xds_parameters.html

Example XSCALE.INP

```
!===== EXAMPLE 3: specific reindexing of input data sets
!
!      Use of specific reindexing of input data sets for resolving
!      indexing ambiguities in the scaled output data set. This
!      happens if the crystal's space group symmetry is lower than
!      its lattice symmetry.
!
RESOLUTION_SHELLS= 100 10 6 4 3 2 1.9
SPACE_GROUP_NUMBER=78
UNIT_CELL_CONSTANTS=57.39 57.39 106.9    90 90 90
OUTPUT_FILE=scaf8_all_merged.hkl
MERGE=TRUE FRIEDEL'S_LAW=FALSE
STRICT_ABSORPTION_CORRECTION=TRUE
INPUT_FILE= ../xds-1_2/XDS_ASCII.HKL
REIDX_ISET= -1  0  0  0  0  1  0  0  0  0 -1  0
INPUT_FILE= ../xds-2_1/XDS_ASCII.HKL
INPUT_FILE= ../xds-3_1/XDS_ASCII.HKL
INPUT_FILE= ../xds-1_4/XDS_ASCII.HKL
INPUT_FILE= *../xds-5_1/XDS_ASCII.HKL
```

Bold keyword/parameter pairs are required. Complete documentation at
xds.mr.mpg.de/html_doc/xscale_parameters.html


Example XDSCONV.INP

```
! UNIT_CELL_CONSTANTS= 10 20 30 90 90 90
! SPACE_GROUP_NUMBER= 96
! GENERATE_FRACTION_OF_TEST_REFLECTIONS=0.05

INPUT_FILE=XDS_ASCII.HKL
OUTPUT_FILE=temp.hkl CCP4_I+F ! or CCP4_I or CCP4_F or SHELX or CNS
FRIEDEL'S_LAW=FALSE ! store anomalous signal in output file even if weak
```

Bold keyword/parameter pairs are required. Complete documentation at xds.mr.mpg.de/html_doc/xdsconv_parameters.html

XDSCC12: calculates $\Delta CC_{1/2,i} = CC_{1/2,with_i} - CC_{1/2,without_i}$



Examples of single data sets (plots from *XDSGUI*)

three resolution ranges (blue=low green=medium red=high) - i refers to batches of width 1°

- find bad frame ranges
- radiation damage

XDSCC12* for selection of data sets for *XSCALE

For important data, do not rely on automatic data processing

There are no clear-cut criteria for the quality of data set (“Table 1” is almost useless).
Synchrotrons typically run multiple pipelines, and the user has to choose ...

What can go wrong in automatic data processing?

- does not handle radiation damage (discard frames towards the end of the data set)
- does not handle shadowed areas of detector
- does not handle indexing problems (multiple lattices, ice, ...) flexibly
- does not optimize processing

Automatic processing with GlobalPhasing's *autoPROC* (and *xia2* for DLS data) is rather reliable for good data sets! **But difficult data sets typically benefit from human insight.**

Manual processing with *XDSGUI*

- problems in phasing and refinement
often due to bad / wrong data processing
- visually inspect frames; mask shadows
- optimize parameters, frame range,
resolution cutoff ..
- presents tables as plots
- interfaces to *XDS* through its files, e.g. **XDS.INP**
- user – extensible / modifiable commands

Manual checks are easily performed with *XDSGUI*

- before running *XDS*: check frames (*adxv*, *Albula*, *XDSGUI* ...)
- **Frame** tab: FRAME.cbf for match obs/calc, and masked areas
- **IDXREF** tab: for ice rings / indexed vs un-indexed reflections
- look at indexing in reciprocal space: **tools** tab/Further analyses/spot2pdb-coot command; see [XDSwiki:XDSGUI#tools](#)
- **statistics** tab: *XDSCC12* and *XDSSTAT*

Optimize data processing

- XDSwiki:Optimisation#Re-INTEGRATEing_with_the_correct_spacegroup.2C_refined_geometry_and_fine-slicing_of_profiles

XDSwiki:Optimisation#using_the_refined_values_for_beam_divergence_and_mosaicity_for_re-integration

- **tools** tab: “Saving and comparing good results” and “Optimizing data quality”. After changing parameters, run “JOB=DEFPIX INTEGRATE CORRECT”, compare and save if better/restore old if worse.
- consider use of *StarAniso* if anisotropy
... may make the difference between interpretable and non-interpretable map, substructure solution, ...

Thank you!