


Multi-crystal approaches


CCP4 Workshop

30.11.2021



Selina Storm



Project Manager for EMBL@PETRA IV





Outline

- motivation
- multi-crystal approaches and serial crystallography
- data collection of multiple crystals - things to consider
- processing data of multiple crystals
 - an introduction to dials.multiplex
 - *in situ* data collection on the main protease (MPro) of COVID-19


Why would one want to use multiple crystals?



- radiation damage
 - especially an issue with small crystals
- to achieve
 - higher resolution and better anomalous signal
 - higher completeness (limitations due to sample mount)
- requirement for time-resolved (serial) experiments





Multiple crystal approaches



Choose your weapon – what are your crystals like?

Adapt your experimental approach accordingly to the aim/situation!


- lots of big, badly diffracting crystals:
collect small sweeps in loops and try to get higher resolution
- brittle crystals which don't like to be taken out of the plate:
grow them in plates which can be mounted on the beamline
(VMXi, I24@Diamond, P14@PETRA III, X06DA@SLS...)





What is the aim of your experiment?

- low-dose / time-resolved structures: serial crystallography (! 'still' images !)

Fixed targets




LCP extruder



Challenges when using multiple crystals

- strictly speaking, each dataset is an experiment in itself
 - variation of the X-ray beam, temperature, humidity, solvent can play a role
- there are no two crystals alike
 - different unit cell parameters
 - different mosaicity
 - different resolution
- preferred orientation (rods, plates)




Things to consider for multi-crystal data collection

- take extra care to keep things as constant as possible (crystallization conditions, temperatures, beamsizes etc)
- vary the starting angle, especially when you have preferred orientation
- key parameter : completeness!
 - ideally, data are processed on the fly
- name your files consistently



Things to consider for multi-crystal analysis


- all data sets need to be processed with the same space group (can be tricky for very small sweeps)
- unit cells need to be quite similar
- data will be scaled together



Programs for data processing – rotational data


- XSCALE (for data processed with XDS)
- BLEND (for data processed with MOSFLM, XDS or DIALS)
- `dials.multiplex` (for data processed with DIALS)
 - generates an html file with an excellent overview
 - autoprocessing on Diamond beamlines to give feedback on-the/fly

Dials.multiplex: how does it work?





Example: structure of M^{pro} binding ligands

- M^{pro}: main protease of SARS-CoV-2
- key enzyme of coronaviruses having a central role in mediating viral replication and transcription
- active site empty and solvent accessible - perfect for fragment screening
- more than 2000 fragments screened, including Moonshot compounds




Example: RT structure of M^{pro} binding ligands

- all initial screens performed on I04-1 at 100 K - are the structures the same at room temperature?
- *in situ* room temperature data collection in plates
- a few drops with the same compound
- strong radiation damage (20 degree)
- preferred orientation, SG C121
 - varying starting angles




Processing on-the-fly

- to answer the question: do we have enough data to get high completeness and sufficient resolution to see the ligand?
- 'force' the correct space group for auto-processing by entering it into ISPyB
- raw data all need to go into the same folder
- two multiplex auto-processing runs in ISPyB:



Processing on-the-fly

- two ISPyB outputs:
 - one processes with the user-defined space group
 - the other determines the space group itself using `dials.cosym` and `dials.symmetry`
- typically 10 data sets needed for a complete data set, more used to increase resolution
- automated molecular replacement with DIMPLE possible



Processing via command line

- required: path to integrated dials results

```
xia2.multiplex $path_to_dials_results.{refl,expt} options
```

- common options:

- `symmetry.space_group=C2`
- `resolution.d_min=2.7`
- `filtering.method=deltacchalf`
`filtering.deltacchalf.stdcutoff=1.5`
- `dose=1,100`

Understanding the results: the multiplex html report

2 3dii ! xia2 dials ! 19x xia2.multiplex ! 19x xia2.multiplex autoPROC+

19 Data Sets Plots Archive Logs & Files

Auto Processing Attachments: 19x xia2.multiplex

Attachments

File	Type	
xia2.multiplex.html	Log	Download View

Anom Multipliciti

6.6
6.4
6.6

Summary

All data

Detailed statistics for dataset All data

Overall

	Overall	Low resolution	High resolution
Resolution (Å)	56.26 - 2.13	56.28 - 5.78	2.17 - 2.13
Observations	105563	5420	5485
Unique reflections	15537	810	801
Multiplicity	6.8	6.7	6.8
Completeness	99.69%	98.78%	100.00%
Mean I/σ(I)	9.6	38.9	0.7
R _{merge}	0.120	0.057	2.168
R _{meas}	0.130	0.062	2.350
R _{pim}	0.048	0.023	0.884
CC _{1/2}	0.996	0.995	0.268

Xia2.multiplex.html

xia2.multiplex repo

Summary

All data

Cosym Analysis

Cosym plots

Cluster comparison

Data sets

Unit cell analysis

Orientation analysis

Delta CC_{1/2} analysis

Intensity clustering

Summary

All data

Detailed statistics for dataset All data

Overall

	Overall	Low resolution	High resolution
Resolution (Å)	56.26 - 2.13	56.28 - 5.78	2.17 - 2.13

Resolution shells

Xtriage

Analysis plots

Analysis by resolution

Analysis by image number

Miscellaneous

Xtriage (all data)

Xtriage

! 1 alert

! The resolution of the data may be lower than the given resolution.

✓ 8 checks passed

Zwart, P. H., Grosse-Kunstleve, R. W. & Adams, P. D. (2005). *CCP4 News*. **43**, contribution 7.

Xtriage

✗ 1 serious warning

✗ The merging statistics indicate that the data may be assigned to the wrong space group.

! 1 alert


! One or more twin operators show a significant twin fraction but since the intensity statistics do not indicate twinning, you may have an NCS rotation axis parallel to a crystallographic axis.

✓ 8 checks passed


Zwart, P. H., Grosse-Kunstleve, R. W. & Adams, P. D. (2005). *CCP4 News*. **43**, contribution 7.

Analysis by resolution (all data)

$\langle I/\sigma(I) \rangle$ vs resolution




$CC_{1/2}$ vs resolution




More analysis by resolution (all data)


Completeness vs resolution



Wilson intensity plot




Analysis by image number - looking for radiation damage



- Scale and Rmerge as well as $I/\sigma(I)$ vs batch plots help to decide which data sets or frames to exclude because of radiation damage

```
xia2.multiplex  
$path_to_dials_results.{refl,expt}  
dose=1,50
```

(to include the first 50 frames per data set)




```
xia2.multiplex  
$path_to_dials_results.{refl,expt}  
exclude_images=7:600:650
```


(to exclude the last 50 images from the 7th data set)

Unit cell analysis (summary)

Distribution of unit cell parameters




Histogram of unit cell parameters



Unit cell analysis (summary)


Unit cell clustering



- can help to find data sets which should be excluded due to non-isomorphism

$\Delta CC_{1/2}$ outlier analysis

Histogram of Delta $CC_{1/2}$




$$\Delta CC_{1/2-i} = CC_{1/2-i} - CC_{1/2\text{-overall}}$$


- excluding the worst negative $\Delta CC_{1/2}$ can improve the overall result

Orientation analysis (in summary)



Stereographic projection (hkl=100)



Stereographic projection (hkl=001)



Example of something going wrong



General advice for processing multiple data sets

- try to process data on-the-fly and keep an eye on completeness and resolution **as a first step!**
- check the analysis carefully for radiation damage and non-isomorphism and identify data sets to cut or exclude
- before you exclude data, check in the orientation analysis whether this data set is critical for completeness
- rerun xia2.multiplex manually and look for improvements

Take-home messages



- plan your experiment and choose your weapon
- collect your data carefully (parameters, path names)
- use auto-processing (xia2.multiplex), but use it wisely

Acknowledgements

Robin Owen
Danny Axford
Sam Horrell

Richard Gildea
James Beilsten-Edmands
Nick Devenish
Graeme Winter

Adam Crawshaw



Thank you for your attention!