

Resolution, data quality and all that ...

Clemens Vornrhein
Global Phasing Ltd.

DLS-CCP4 Workshop 2021

Too much focus on single number (resolution) to describe model quality?

intitle:resolution ccp4bb site:www.mail-archive.com

All News Videos Shopping Maps More Settings

About **564** results (0.34 seconds)

intitle:diffraction ccp4bb site:www.mail-archive.com

All Videos Shopping News Images More Settings

About **283** results (0.39 seconds)

intitle:quality ccp4bb site:www.mail-archive.com

All Maps Shopping News Videos More Settings

About **58** results (0.31 seconds)

x 2

x 10


x 5

Resolution = ability to resolve detail

1.0 Å

2CNQ (1.00 Å)
=13.5

B(A244)=7.3




1.5 Å

6N19 (1.50 Å)
=27.0

B(A522)=17.6


1.4 Å



2.0 Å

6J66 (1.95 Å)
=28.8

B(A181)=27.5




2.5 Å

6NT2 (2.48 Å)
=28.4

B(B216)=17.9

2.5 Å




3.0 Å

6JJ9 (3.00 Å)
=56.0

B(A201)=31.5


3.3 Å



3.5 Å


6DC5 (3.50 Å)
=72.2

B(I35)=31.0




2mFo-DFc maps after BUSTER refinement at 1.0 rms

Resolution, diffraction limit and binning



- ❑ main author: Ian Tickle (Global Phasing)
- ❑ **anisotropic** cut-off of merged intensity data
 - avoids problems of
 - isotropic resolution cut-offs
 - throwing away good data (in one direction)
 - including poor data (in another direction)
- ❑ Bayesian estimation of structure factor amplitudes
 - more rigorous treatment of prior and anomalous data
 - for more details see:
https://staraniso.globalphasing.org/staraniso_about.html
- ❑ (optional) anisotropic correction of data
- ❑ Part of autoPROC ('process', 'aP_scale')



merged or
unmerged data
(MTZ, XDS)

Anisotropy is not necessarily ellipsoidal!

Isotropic vs anisotropic

GLOSSARY of terminology related to anisotropy

ABOUT ANISOTROPY

Anisotropy is the property of being directionally dependent, implying a variation of certain physical properties as a function of direction. It is the opposite of **isotropy**, that is, possessing exactly the same properties in whatever direction one looks. Here we are mainly concerned with **two** potentially anisotropic properties of the measured X-ray diffraction data: the **diffraction limit** (also known as the '**resolution limit**'), and the **attenuation of the diffracted intensity** due to thermal motion and/or disorder.

The diffraction limit

Beyond the diffraction limit of the sample, the merged intensities of unique reflections fall, on average, below the level of statistical significance and should therefore be considered to be non-informative (inclusion of such data in further procedures, scaling in particular, may in fact be positively detrimental).


The statistical significance of the intensity data after merging is commonly expressed as either the value of the mean signal/noise ratio $I/\sigma(I)$, or the half-dataset correlation coefficient $CC_{1/2}$ [1], computed in spherical shells of reciprocal space. The statistical significance, in whatever form one chooses to express it, is used as a criterion for applying a **resolution cut-off** to the reflection data that should match the diffraction limit as closely as possible. Clearly if the applied resolution cut-off falls inside the diffraction limit, then useful data containing information about the structure will be lost; if it falls outside, then uninformative data will be kept.

Since the diffraction limit may be anisotropic, in STARANISO the locally-averaged value of $I/\sigma(I)$, rather than its spherical average, is used to define an anisotropic resolution cut-off - that is, a cut-off varying with direction in reciprocal space. The boundary between statistically significant and insignificant data therefore defines a surface, the **anisotropic resolution cut-off surface** (see diagram on right: click image to enlarge).

No assumption is made about the shape or orientation of the anisotropic resolution cut-off surface, other than that it has at least the symmetry of the Laue class. This means that in orthorhombic and higher symmetry space groups the principal axes of the surface shape must coincide with the crystal axes. In monoclinic space groups only one of the principal axes of the surface shape is constrained to lie along a crystal axis, namely **b**, but there is no symmetry constraint on the second principal axis: it can lie anywhere in the **ac** plane since the choices of the **a** and **c** axes are completely arbitrary. Similarly, in triclinic space groups the choices of all three crystal axes are arbitrary so there are no symmetry constraints whatsoever on the directions of the principal axes of the surface shape.

Note that it is **not** assumed that the shape of the anisotropic cut-off surface is an ellipsoid (the contours in the figure below are clearly not ellipsoidal). The principal cause of distortion of the diffraction limit surface away from a perfect ellipsoid is the variation of the reflection redundancy through its effect on the $\sigma(I)$ values of the merged data (*i.e.* inversely proportional to the square root of the redundancy), and therefore on the mean $I/\sigma(I)$. The variation of redundancy in reciprocal space is largely arbitrary, being a consequence of both the experimental setup and the data collection strategy (whether a standard protocol or a deliberate user choice).


However, for the purposes of estimating the resolution limit as a function of direction an ellipsoid is fitted to the anisotropic cut-off surface. Measured data with d^* greater than the largest semi-axis length of the fitted ellipsoid are considered to be outside the anisotropic cut-off surface. Then all measured intensities (excluding any systematic absences) that fall inside this anisotropic cut-off surface are considered to be '**observed**', that is they are included in the reflection statistics and are written to the reflection output file, while all those that fall outside are '**unobserved**' and are not used further. 'Unmeasured data (except systematic absences)' that fall inside the fitted ellipsoid (for example in a cusp) are deemed to be '**observable**' (*i.e.* they had the potential to be observed if only they had been measured), whereas those unmeasured reflections falling outside the fitted ellipsoid are deemed to be '**unobservable**' (*i.e.* not worth measuring).




autoPROC: STARANISO plots

rec. unit cell

Unmeasured and (expected to be) unobservable
Unmeasured and (expected to be) observable
Measured and unobserved

Diffraction limits
observableDiffraction limits
recorded

3D WebGL visualisation in STARANISO


<https://staraniso.globalphasing.org/cgi-bin/PDBpeep.cgi>
<https://staraniso.globalphasing.org/>
local $\langle I/\sigma I \rangle$

4YUO P_{2,2,2}_{1,1,1}
 $a^* \sim 7^\circ$
 single sweep


beam centre in
 module gap

merged data


unmerged data



redundancy

local $\langle I/\sigma I \rangle$ 


4ZC9 P_{2,2,2}_{1,1,1}
 $b^* \sim 14^\circ$
 2 sweeps
 (distance)



- **The following results were communicated to us by Neil Paterson at Diamond Light Source.**

The panels show the maps after density modification (CCP4 DM program) using (left) the standard protocol with an isotropic resolution cut-off, and (right) the STARANISO protocol with anisotropic resolution cut-off and correction. Compare in particular the density in the more mobile loop regions at the ends of the helices, and also the flatness of the solvent region (right-hand half of each map).

Neil comments: "STARANISO has been key to solving our structure. We've got a severely anisotropic case that goes from around 7 Å along a^* to 3.0-3.5 Å along b^*/c^* , with lots of artifacts in the maps. We can find the substructure using Shel-X without correction but the maps are pretty bad; once corrected it works a treat (77% solvent!)."



Evolution of diffraction data metrics - “Table 1”

 $R_{\text{merge}} (R_{\text{sym}})$

Completeness

Redundancy (multiplicity)

 $\langle I/\sigma \rangle$ R_{meas} R_{pim} $CC_{1/2}$

As function of resolution in spherical shells (equal spacing or equal volume): overall and outer shell

Table I. X-ray Data Reduction and Crystallographic Refinement Statistics

(A) X-ray data reduction statistics	
Space group	$P2_12_12_1$
Unit cell dimensions	
a, b, c	33.5 Å, 68.2 Å, 137.6 Å
Resolution	68.8 Å–2.25 Å
(last shell)	(2.26 Å–2.25 Å)
Total measurements (last shell)	109,559 (1,168)
Number of unique reflections	15,637
(last shell)	(163)
Wavelength	1 Å
R_{merge}^a (last shell)	0.065 (0.592)
$I/\sigma(I)$ (last shell)	16.4 (3.3)
Completeness (last shell)	0.997 (1)
Redundancy (last shell)	7 (7.2)

R_{merge} : Arndt et al (1958)

R_{meas} : Diederichs & Karplus (1997), Weiss & Hilgenfeld (1997)

R_{pim} : Weiss (2001)

$CC_{1/2}$: Karplus & Diederichs (2012)

❑ Equally spaced in $(d^*)^2$

- Number of measurements used for computing bin averages increases with resolution
- Low-resolution statistics can be unreliable (too few measurements)
- High-resolution statistics not finely enough sampled?
- Used by XDS, AIMLESS

SUBSET OF INTENSITY DATA WITH SIGNAL/NOISE ≥ -3.0

RESOLUTION LIMIT	NUMBER OF REFLECTIONS			COMPLETENESS OF DATA
	OBSERVED	UNIQUE	POSSIBLE	
3.13	75353	14441	14448	100.0%
2.22	136754	26193	26370	99.3%
1.81	178487	34010	34139	99.6%
1.57	207171	40372	40378	100.0%
1.40	214727	45702	45725	99.9%
1.28	183266	50533	50658	99.8%
1.18	142566	54470	55100	98.9%
1.11	57895	43932	59148	74.3%
1.04	15808	15171	63013	24.1%
total	1212027	324824	388979	83.5%

❑ Equal volume (to have same number of reciprocal lattice points in each bin)

- Allows different sampling (coarse or fine) depending on requirements
- Remember ice-rings (3.89 3.66 3.44 2.67 2.25 2.07 1.95 1.91 1.88 1.72 1.52 1.47 1.44 1.37 1.36 1.30 1.27 1.26 1.22 1.17 1.12)!
- Adequate for isotropic data that is complete with homogeneous multiplicity

❑ Equal number of (actual) observations

- Can be seen as generalisation of idea behind “equal volume” binning
- Automatically self-adjusting for anisotropic (STARANISO) and incomplete data (serial crystallography, LCP, in-plate, ...)

Binning: decision making, data description

Binning is boring - but an important aspect here:

☐ Decision making:

- Which images/datasets/reflections to include
- Which reflections to include: isotropic/anisotropic diffraction limit

☐ Comparisons:

- Between programs and pipelines
- Between different processing options

☐ Basis of decision making


- Raw data (binned statistics)
- Smoothed (spline, Bezier, ...)
- Ice-ring resolution ranges included/excluded/smoothed?

What is my high-resolution limit? ... depending on binning method


d^*^2 (20 bins)



equal-volume





equal-Nobs



$I/sigI > 1$

looking from the low-resolution end






How do we look at data/statistics: from low- or high-resolution end?

d^*^2 (20 bins)




equal-volume

equal-Nobs




CC(1/2)>0.3

looking from the high-resolution end



automatic processing at synchrotrons



Some generalisation required to provide similar/identical information for different pipelines.

Be aware of different binning (ie don't over-interpret small differences)

OSC 17-05-2018 16:30:00
 /data/visitor/mx415/ld30a2/20180517/RAW_DATA/4

Pipeline	SpaceGroup	a (Å)	b (Å)	c (Å)	α (°)	β (°)	γ (°)	Shell	Resolution (Å)	Multiplicity	Completeness %	<I/Sigma>	Rmeas	Rmerge	Rpim	cc(1/2)	ccAno	sigAno	ISA	Download
BEST autoPROC	123	79.0	79.0	79.0	90.0	90.0	90.0	Overall Inner Outer	55.9-2.3 55.9-6.3 2.4-2.3	5.2 5.3 5.4	99.1 94.7 100.0	16.0 36.3 2.8	8.3 2.6 47.7		100 100 90					
autoPROC_staranisotropy	123	79.0	79.0	79.0	90.0	90.0	90.0	Overall Inner Outer	55.9-1.8 55.9-5.1 1.9-1.8	5.1 5.0 5.3	92.0 97.0 49.0	9.6 34.7 -0.2	10.5 3.1 -607.4		100 100 40					
BEST autoPROC	123	79.0	79.0	79.0	90.0	90.0	90.0	Overall Inner Outer	55.9-2.3 55.9-6.3 2.4-2.3	5.2 5.3 5.5	98.8 94.7 99.5	16.1 36.7 2.7	9.1 2.6 51.8		100 100 90					
autoPROC_staranisotropy	123	79.0	79.0	79.0	90.0	90.0	90.0	Overall Inner Outer	55.9-1.8 55.9-5.2 1.9-1.8	5.1 5.0 5.4	91.1 97.0 48.9	10.0 35.0 -0.6	11.6 3.2 -168.9		100 100 60					

also: ALBA (ES), APS (US), MaxIV (S), PETRA (D), SLS (CH), SSRF, ...

Diamond (DLS), UK

ESRF, F

Rupp, B., 2018. Against Method: Table 1—Cui Bono?. Structure, 26(7), pp.919-923.

Structure

Perspective

Against Method: Table 1—Cui Bono?


Bernhard Rupp^{1,2,*}

¹k.-k.Hofkristallamt, San Diego, CA 92084, USA

²Division of Genetic Epidemiology, Medical University Innsbruck, Schöpfstraße 41, Innsbruck, Tyrol 6020, Austria

*Correspondence: bernhard.rupp@i-med.ac.at

<https://doi.org/10.1016/j.str.2018.04.013>



Diffraction data quality metrics of **PDB** entries harvested from archived mmCIF files - from depositions (20211028 PDB release, aB_cif2table1 from **BUSTER** refinement package, **Global Phasing Ltd.**)

PDB identifier (4-letter code):

Go

- by PDB id
- by data collection year
- by Synchrotron Beamline
- by Data scaling software
- by Refinement software

You can get to a particular entry (e.g. 5UHO) via the URL

<https://staraniso.globalphasing.org/table1/uh/5uho.html>. This resource is automatically updated on a weekly basis. For any questions, please contact us at **proc-devel@globalphasing.com**.

Last modification: Thu Oct 28 05:38:42 CEST 2021

General information		
Spacegroup name	P 32 2 1	
Unit cell parameters	47.914 47.914 200.488 90.0 90.0 120.0	
Wavelength	0.97932 Å	
Data quality metrics		
	Overall	OuterShell
Low resolution limit [Å]	40.630	1.940
High resolution limit [Å]	1.900	1.900
Rmerge	0.050	0.433
Rmeas	0.062	0.540
Rpim	0.036	0.318
Total number of observations	-	
Total number unique	22045	1386
$\langle I/\sigma(I) \rangle$	14.30	3.10
Completeness [%]	99.9	99.9
Multiplicity	4.8	4.8
CC(1/2)	0.997	0.913

<https://staraniso.globalphasing.org/table1/>

“Table 1” - 5E74 : isotropic (standard)

Spacegroup name C2221
 Unit cell parameters 80.694 96.593 57.922 90.0 90.0 90.0
 Wavelength 0.97942 Å

“Table 1” from autoPROC

Number of active ice-rings within this resolution range = 7

Criteria used in determination of resolution cut:

 Rpim <= 0.6000
 I/sigI >= 2.00
 CC(1/2) >= 0.3000

assuming isotropic data

additional items:

- active ice-rings
- criteria for cut-off decision
- ... others?

	Overall	InnerShell	OuterShell
Low resolution limit	42.302	42.302	1.686
High resolution limit	1.657	4.497	1.657
Rmerge (all I+ & I-)	0.035	0.020	0.819
Rmerge (within I+/I-)	0.035	0.019	0.814
Rmeas (all I+ & I-)	0.038	0.021	0.873
Rmeas (within I+/I-)	0.039	0.022	0.927
Rpim (all I+ & I-)	0.013	0.007	0.297
Rpim (within I+/I-)	0.018	0.010	0.438
Total number of observations	204143	12396	11170
Total number unique	23978	1462	1325
Mean(I)/sd(I)	25.2	74.7	2.3
Completeness	88.0	99.9	99.9
Multiplicity	8.5	8.5	8.4
CC(1/2)	1.000	1.000	0.931
Anomalous completeness	87.4	99.9	99.8
Anomalous multiplicity	4.5	4.9	4.4
CC(ano)	-0.245	-0.208	-0.025
DANO /sd(DANO)	0.651	0.689	0.572

"Table 1" - 5E74 : anisotropic (STARANISO)

Spacegroup name C2221
 Unit cell parameters 80.694 96.593 57.922 90.0 90.0 90.0
 Wavelength 0.97942 Å

Diffraction limits & eigenvectors of ellipsoid fitted to diffraction cut-off surface:

2.186	1.0000	0.0000	0.0000	\hat{a}^*
1.442	0.0000	1.0000	0.0000	\hat{b}^*
1.379	0.0000	0.0000	1.0000	\hat{c}^*

Number of active ice-rings within this resolution range = 7

Criteria used in determination of diffraction limits:

local(I/sigI) >= 1.20


	Overall	InnerShell	OuterShell
Low resolution limit	42.302	42.302	1.504
High resolution limit	1.386	4.824	1.386
Rmerge (all I+ & I-)	0.035	0.019	1.146
Rmerge (within I+/I-)	0.034	0.019	1.074
Rmeas (all I+ & I-)	0.037	0.020	1.230
Rmeas (within I+/I-)	0.038	0.021	1.234
Rpim (all I+ & I-)	0.012	0.007	0.435
Rpim (within I+/I-)	0.018	0.009	0.596
Total number of observations	202888	10343	9127
Total number unique	24057	1202	1203
Mean(I)/sd(I)	25.5	75.5	1.4
Completeness (spherical)	52.0	100.0	12.1
Completeness (ellipsoidal)	84.1	100.0	72.9
Multiplicity	8.4	8.6	7.6
CC(1/2)	1.000	1.000	0.645
Anomalous completeness (spherical)	50.8	99.9	11.5
Anomalous completeness (ellipsoidal)	83.3	99.9	72.3
Anomalous multiplicity	4.5	5.0	4.0
CC(ano)	-0.213	-0.229	-0.015
DANO /sd(DANO)	0.656	0.679	0.629

Diffraction
limits and
directions

Completeness:


#(HKL with value) /
#(possible HKL)

#observed /
#observable




Spots indexed vs image number per unique orientation (“iterative indexing”)

4XGU




5J1S




How to report this in 'Table 1'?


5SZB




5OPE



STARANISO within autoPROC – 2D plots



ice



cusp + ice

gaps


cusp

Unindexed spots on detector surface ice or powder rings


How do we report and record this:

- to the user (to be aware of e.g. cryo-cooling problems)
- upon deposition (to alert user of structure to potential contamination of data)


5E73
1.71Å




5HQG
2.00Å





5USC
2.20Å



6CKP
1.15Å




Summary HTML (stages, decision making)



ISPyB XML file

staraniso_alldata-unique.mtz
staraniso_alldata-unique.table1
Data_1_autoPROC_STARANISO_all.cif

PDF report



- ❑ Is there such a thing as “the resolution”?
 - ❑ We should probably start talking about diffraction limit(s).

 - ❑ Describing diffraction data via a single number makes life easy ... but doesn't really do the variation between different datasets justice: multiple lattices, ice, detector gaps, cusp, split spots etc.

 - ❑ Thinking about binning is boring - but important when assuming isotropic diffraction and making decisions.

 - ❑ The STARANISO approach (local $\langle I/\sigma(I) \rangle$ to define anisotropic cut-off surface) excludes noise and includes all significant data while avoiding (isotropic) binning headache.
 - Completeness: fraction of data actually observed relative to potentially observable data (i.e. significant, not just Miller indices).
-

Resolution, data quality and all that ...

Clemens Vonnrhein
Global Phasing Ltd.

DLS-CCP4 Workshop 2021

C. Vonnrhein, Global Phasing Ltd
DLS-CCP4 2021

Too much focus on single number (resolution) to describe model quality?




C. Vornrhein, Global Phasing Ltd
DLS-CCP4 2021

Especially resolution seems to preoccupy users a lot: it is mentioned 10 times more often in the subject line of messages to the CCP4BB than e.g. the term “quality”.

Why is that?


Because we traditionally equate data resolution with model quality: the higher the resolution, the better the model.

Resolution = ability to resolve detail




G. Vonnheim, Global Phasing Ltd
DLS-CCP4 2021

Resolution, diffraction limit and binning



STARANISO – <http://staraniso.globalphasing.org/>

- ☐ main author: Ian Tickle (Global Phasing)
- ☐ **anisotropic** cut-off of merged intensity data
 - avoids problems of
 - isotropic resolution cut-offs
 - throwing away good data (in one direction)
 - including poor data (in another direction)
- ☐ Bayesian estimation of structure factor amplitudes
 - more rigorous treatment of prior and anomalous data
 - for more details see:
https://staraniso.globalphasing.org/staraniso_about.html
- ☐ (optional) anisotropic correction of data
- ☐ Part of autoPROC ('process', 'aP_scale')
 - user is given classical, isotropic (TRUNCATE) output
 - but also data after **anisotropic** analysis (STARANISO)



merged or
unmerged data
(MTZ, XDS)

Anisotropy is not necessarily ellipsoidal!

Isotropic vs anisotropic

GLOSSARY of terminology related to anisotropy

ABOUT ANISOTROPY

Anisotropy is the property of being directionally dependent, implying a variation in certain physical properties as a function of direction. It is the opposite of isotropy, that is, possessing exactly the same properties in whatever direction you look. That is, an anisotropic material has directional physical properties, while isotropic materials have the same properties in all directions.

The diffraction limit

Figure 1: Diffraction limit. The diffraction limit is the smallest resolution of a system, determined by the wavelength of the light used. It is the limit beyond which the system cannot resolve details smaller than the wavelength of the light used.


The resolution of a system is the smallest distance between two points that can be distinguished. It is determined by the wavelength of the light used and the numerical aperture of the system. The resolution of a system is the smallest distance between two points that can be distinguished.

There are different ways to define resolution. The most common is the Rayleigh criterion, which states that two points are resolved if the central maximum of one diffraction pattern falls on the first minimum of the other. Other criteria include the Sparrow criterion and the Dawes criterion.

In anisotropic systems, the resolution is not the same in all directions. This is because the system's properties vary with direction. For example, in a crystal, the resolution is higher along the crystallographic axes than in other directions. This is due to the anisotropic nature of the crystal's structure.


There are several ways to define resolution. The most common is the Rayleigh criterion, which states that two points are resolved if the central maximum of one diffraction pattern falls on the first minimum of the other. Other criteria include the Sparrow criterion and the Dawes criterion.

There are several ways to define resolution. The most common is the Rayleigh criterion, which states that two points are resolved if the central maximum of one diffraction pattern falls on the first minimum of the other. Other criteria include the Sparrow criterion and the Dawes criterion.



C. Vonrhein, Global Phasing Ltd
DLS-CCP4 2021

autoPROC: STARANISO plots




C. Vornheim, Global Phasing Ltd
DLS-CCP4 2021

3D WebGL visualisation in STARANISO

<https://staraniso.globalphasing.org/cgi-bin/PDBpeep.cgi>

<https://staraniso.globalphasing.org/>




C. Vornrhein, Global Phasing Ltd
DLS-CCP4 2021

• The following results were communicated to us by Neil Paterson at Diamond Light Source.

The panels show the maps after density modification (CCP4 DM program) using (left) the standard protocol with an isotropic resolution cut-off, and (right) the STARANISO protocol with anisotropic resolution cut-off and correction. Compare in particular the density in the more mobile loop regions at the ends of the helices, and also the flatness of the solvent region (right-hand half of each map).

Neil comments: 'STARANISO has been key to solving our structure. (We've got a severely anisotropic case that goes from around 3.4 Å along a* to 3.0-3.5 Å along b* (c*, with lots of artifacts in the maps. We can find the substructure using Shel-X without correction but the maps are pretty bad. Once corrected it works a treat (77% solvent)).'



G. Vornrhein, Global Phasing Ltd
DLS-CCP4 2021

Evolution of diffraction data metrics - "Table 1"

$R_{\text{merge}} (R_{\text{sym}})$

Completeness

Redundancy (multiplicity)

$\langle I/\sigma \rangle$

R_{meas}

R_{pim}

$CC_{1/2}$

As function of resolution in spherical shells (equal spacing or equal volume): overall and outer shell

Table 1. X-ray Data Reduction and Crystallographic Refinement Statistics

(A) X-ray data reduction statistics	
Space group	$P2_12_12_1$
Unit cell dimensions	
a, b, c	33.5 Å, 68.2 Å, 137.6 Å
Resolution	68.8 Å–2.25 Å
(last shell)	(2.26 Å–2.25 Å)
Total measurements (last shell)	109,558 (1,168)
Number of unique reflections	15,637
(last shell)	(163)
Wavelength	1 Å
R-merge* (last shell)	0.065 (0.592)
I/σ(I) (last shell)	16.4 (3.3)
Completeness (last shell)	0.997 (1)
Redundancy (last shell)	7 (7.2)

R_{merge} : Arndt et al (1958)

R_{meas} : Diederichs & Karplus (1997), Weiss & Hilgenfeld (1997)

R_{pim} : Weiss (2001)

$CC_{1/2}$: Karplus & Diederichs (2012)

C. Vornrhein, Global Phasing Ltd
DLS-CCP4 2021

The metrics for describing data quality have evolved over the last few years

New variants of R-values have appeared and the use of correlation coefficients is now standard.

Our picture of how to best describe data has clearly evolved – so maybe it is time to also look at our notion of resolution and completeness?

Binning

□ Equally spaced in $(d^*)^2$

- Number of measurements used for computing bin averages increases with resolution
- Low-resolution statistics can be unreliable (too few measurements)
- High-resolution statistics not finely enough sampled?
- Used by XDS, AIMLESS

Summary of intensity data with global/unique $\sigma = -3.0$

RESOLUTION LIMIT	NUMBER OF OBSERVED	NUMBER OF REFLECTIONS UNIQUE	POSSIBLE	COMPLETENESS OF DATA
2.13	72253	14441	14448	100.0%
2.22	136754	26193	26370	99.3%
1.81	178487	34010	34139	99.6%
1.57	207171	40372	40378	100.0%
1.40	214727	42702	42725	99.9%
1.28	182206	35533	35638	99.8%
1.18	142546	24470	24530	98.9%
1.11	57890	42932	39148	74.2%
1.04	15808	15171	63013	24.1%
LOLO	1212027	324824	388979	83.5%

□ Equal volume (to have same number of reciprocal lattice points in each bin)

- Allows different sampling (coarse or fine) depending on requirements
- Remember ice-rings (3.89 3.66 3.44 2.67 2.25 2.07 1.95 1.91 1.88 1.72 1.52 1.47 1.44 1.37 1.36 1.30 1.27 1.26 1.22 1.17 1.12)!
- Adequate for isotropic data that is complete with homogeneous multiplicity

□ Equal number of (actual) observations


- Can be seen as generalisation of idea behind "equal volume" binning
- Automatically self-adjusting for anisotropic (STARANISO) and incomplete data (serial crystallography, LCP, in-plate, ...)

Binning: decision making, data description


Binning is boring - but an important aspect here:

- ☐ Decision making:
 - Which images/datasets/reflections to include
 - Which reflections to include: isotropic/anisotropic diffraction limit
 - ☐ Comparisons:
 - Between programs and pipelines
 - Between different processing options
 - ☐ Basis of decision making
 - Raw data (binned statistics)
 - Smoothed (spline, Bezier, ...)
 - Ice-ring resolution ranges included/excluded/smoothed?
-


What is my high-resolution limit? ... depending on binning method



C. Vornrhein, Global Phasing Ltd
DLS-CCP4 2021



How do we look at data/statistics: from low- or high-resolution end?



C. Vonrhein, Global Phasing Ltd
DLS-CCP4 2021

Rupp, B., 2018. Against Method: Table 1—Cui Bono?. Structure, 26(7), pp.919-923.

Structure
Perspective

Against Method: Table 1—Cui Bono?

Bernhard Rupp^{1,2,*}

¹UCLA, 621 Charles E. Young Drive South, Los Angeles, CA 90095, USA

²Division of Genetic Epidemiology, Medical University Innsbruck, Schöpfstraße 41, Innsbruck, Tyrol 6020, Austria

*Correspondence: bernhard.rupp@med.ac.at

<https://doi.org/10.1016/j.str.2018.04.013>

Diffraction data quality metrics of PDB entries harvested from archived mmCIF files - from depositions (20211028 PDB release, aB_cif2table1 from BUSTER refinement package, Global Phasing Ltd.)

PDB identifier (4-letter code): Go


- by PDB id
- by data collection year
- by Synchrotron Beamline
- by Data scaling software
- by Refinement software

You can get to a particular entry (e.g. 5UHQ) via the URL <https://staraniso.globalphasing.org/table1/5uho.html>. This resource is automatically updated on a weekly basis. For any questions, please contact us at prc-develop@globalphasing.com.

Last modification: Thu Oct 28 05:38:42 CEST 2021

General information		
Deposition name	PDB	
Deposition date	2021-10-28	
Deposition ID	5UHQ	
Data quality metrics		
	Overall	Overall
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00
Resolution (Å)	4.00	4.00

<https://staraniso.globalphasing.org/table1/>



When it comes to archiving, it is not only the storing of a rich set of meta-data that is important, but also the way of displaying the stored data quality metrics to the user of a particular structural model. No one will have the time to look through long columns of raw numbers. At the same time, an assessment of data quality to judge the quality of the archived model in the context of the experimental data is crucial - especially for users that will have difficulty interpreting the details of an actual electron density map.

In his 2018 paper, Bernhard Rupp argued for a revamping of the traditional “Table 1” and cited the way STARANISO visualises 3D data quality metrics as an alternative. As a way of standardising the already available data quality metrics information of existing PDB entries, we provide a web-service at the URL shown below. This might show a “Table 1”

“Table 1” - 5E74 : isotropic (standard)

Experiment name: 5E74
Unit cell parameters: 88.894 93.894 87.922 98.0 99.0 99.0
Resolution: 0.3042 Å

“Table 1” from autoPROC

Number of active ice-rings within this resolution range = 7

assuming isotropic data

Criteria used in determination of resolution limit:

Spin $\Delta\mu$ 0.6000
C/edge $\Delta\mu$ 2.00
CC1/2 $\Delta\mu$ 0.3000

additional items:

- active ice-rings
- criteria for cut-off decision
- ... others?

	overall	inner-shell	outer-shell
Resolution limit	0.302	0.300	0.306
High resolution limit	0.657	0.693	0.657
Score (all $h \leq 1$)	0.035	0.022	0.019
Score (within $h \leq 1$)	0.035	0.015	0.014
Score (all $h \leq 1$)	0.038	0.021	0.013
Score (within $h \leq 1$)	0.039	0.023	0.027
Spin (all $h \leq 1$)	0.013	0.004	0.002
Spin (within $h \leq 1$)	0.038	0.012	0.008
Total number of observations	25513	12295	11170
Total number unique	23978	1462	1325
R _{meas} ($h \leq 1$)	28.2	74.1	2.3
Completeness	88.0	99.5	99.9
Multiplicity	8.5	8.5	8.4
CC1/2	1.000	1.000	0.931
Anomalous completeness	87.4	99.5	99.8
Anomalous multiplicity	4.5	4.5	4.4
CC _{anom}	+0.245	+0.205	+0.025
SAD/2D/3D	0.451	0.680	0.542

C. Vornheim, Global Phasing Ltd
DLS-CCP4 2021

“Table 1” - 5E74 : anisotropic (STARANISO)

Experiment namec3221Unit cell parameters80.694 95.893 57.922 90.0 90.0 90.0Wavelength0.97942 Å

Diffraction limits & eigenvectors of ellipsoid fitted to diffraction cut-off surface:
2.1861.00000.00000.0000a₁
1.4320.00001.00000.0000b₁
1.3790.00000.00001.0000c₁

Number of active bins/zone within this resolution range = 7

Criteria used in determination of diffraction limits:
Local(I/σ(I)) σ² = 1.00

OverallInnerShellOuterShell

Resolution limit42.10242.3024.894

High resolution limit1.3860.8241.386

Sweep (a1 1=1 1=1)

Score (within 2σ)0.0250.0181.116

Score (all 1=1 1=1)0.0240.0181.074

Score (all 1=1 1=1)0.0270.0201.230

Score (within 2σ)0.0380.0211.234

Score (all 1=1 1=1)0.0120.0040.435

Score (within 2σ)0.0100.0040.596

Total number of observations202000103519127

Total number unique2505712031253

Mean(I)/σ(I)25.525.51.0

Completeness (spherical)92.0100.012.1

Completeness (ellipsoidal)84.1100.072.9

Multiplicity8.48.57.6

CC(1/2)1.0001.0000.645

Completeness:
#(HKL with value) /
#(possible HKL)
#observed /
#observable

Observed data

Observed data

Observed data


Diffraction
limits and
directions

Completeness:


#(HKL with value) /
#(possible HKL)

#observed /
#observable

Spots indexed vs image number per unique orientation ("iterative indexing")




How to report this in "Table 1"?







C. Vornrhein, Global Phasing Ltd
DLS-CCP4 2021

STARANISO within autoPROC – 2D plots



C. Vonrhein, Global Phasing Ltd
DLS-CCP4 2021

Unindexed spots on detector surface ice or powder rings




How do we report and record this:

- to the user (to be aware of e.g. cryo-cooling problems)
- upon deposition (to alert user of structure to potential contamination of data)

C. Vornheim, Global Phasing Ltd
DLS-CCP4 2021


autoPROC reporting

Summary HTML (stages, decision making)



ISPyB XML file

```
stارانiso_alldata-unique.mtz
stارانiso_alldata-unique.table1
Data 1 autoPROC STARANISO all.cif
```

PDF report

Summary

- ☐ Is there such a thing as “the resolution”?
 - ☐ We should probably start talking about diffraction limit(s).

 - ☐ Describing diffraction data via a single number makes life easy ... but doesn't really do the variation between different datasets justice: multiple lattices, ice, detector gaps, cusp, split spots etc.

 - ☐ Thinking about binning is boring - but important when assuming isotropic diffraction and making decisions.

 - ☐ The STARANISO approach (local $\langle I/\sigma(I) \rangle$ to define anisotropic cut-off surface) excludes noise and includes all significant data while avoiding (isotropic) binning headache.
 - Completeness: fraction of data actually observed relative to potentially observable data (i.e. significant, not just Miller indices).
-

C. Vonnrhein, Global Phasing Ltd
DLS-CCP4 2021