

GRAEME WINTER

CCP4 / DIAMOND WORKSHOP 2022

# DATA COLLECTION STRATEGIES FOR PHOTON COUNTING DETECTORS



# OVERVIEW


- Before we start advice: if you remember two things...
- Data collection strategies
- Detectors
- Photon counting detectors
- Diffraction data collection advice
- Tactical suggestions
- Conclusions

IF YOU ARE GOING TO REMEMBER ONE THING...

Ask your beamline scientist / local contact for advice

“What would you suggest?”

AND IF YOU CAN REMEMBER TWO THINGS



DON'T BE DARTH



# DATA COLLECTION STRATEGIES

# GOALS FOR DATA COLLECTION

- Measure all the unique reflections
- Measure enough copies that you can assess how well they agree / determine symmetry etc.
- Maximise the “quality” of the signal you record
- Minimise the radiation damage (quantitative discussion of this will follow in subsequent lectures)







# UNIQUE REFLECTIONS: GEOMETRIC STRATEGY

- Asymmetric unit: whole volume of reciprocal space "divided by" the crystal symmetry - you want all of this ideally
- For high symmetry crystals can get away with a relatively small rotation to get complete data
- For low symmetry data may need much more
- More measurements of reflections also improves the quality of the average




# ASYMMETRIC UNITS






# BLIND REGIONS - ISSUE IN LOW SYMMETRY




# BLIND REGIONS - ISSUE IN LOW SYMMETRY





# BLIND REGIONS - MULTIPLE ORIENTATIONS



# RULES OF THUMB - ASSUMING PIXEL ARRAY DETECTOR

- Unless you have a specific reason not to, collect no less than 360° of data
- If you can, collect with more than one sample orientation w.r.t. beam
- Merging multiple data sets is easy and has ancillary benefits - better coverage of the absorption surface, truly independent measurements of reflections etc.
- Cost of more data with a lower dose is simply in processing time - and how long do you spend getting your samples



# LEGITIMATE REASONS FOR $< 360^\circ$

- Sample mounting precludes full rotation - in situ / many samples on mesh
- I know precisely what I am doing for this fragment screen -  $120 / 180^\circ$  is known to be enough





DETECTORS

# DETECTORS

- Convert X-ray photons into some sort of “intensity value”
- Started off with photographic film (data processing involved chemicals, darkroom, film scanner)
- Most likely to encounter CCD or photon counting pixel array detectors
- No detector is perfect but they are very good - knowing the properties will help you to collect good data



# DETECTORS






# COMMON DETECTOR TYPES

- Integrating:
  - Not worried about photons / second, total photons / pixel problem
  - Read-out noise, dark current
  - Potentially slow (CCD - *not* Jungfrau) - 0.5s / frame read out time
  - Variable gain (counts / photon)
- Photon counting:
  - Limited photons / second rate
  - No read-out noise, dark current
  - Substantially faster - allows shutterless collection
  - Nuances relating to charge sharing etc.



# PHOTON COUNTING DETECTORS





# PHOTON COUNTING DETECTORS


- Direct conversion of photon energy to charge in sensor - derived from sensors from particle physics
- Readout electronics for every pixel (e.g. unlike CCD)
- Enables shutterless data collection as readout is (near-)instant (ms to  $\mu$ s)
- Dectris PILATUS and EIGER most common

# PHOTON COUNTING DETECTORS




Detector: Si sensor material connected  
to readout electronics - large potential  
difference across sensor

# PHOTON COUNTING DETECTORS





# PHOTON COUNTING DETECTORS



X-ray ionizes Si atom


# PHOTON COUNTING DETECTORS



e- cascade

(3500 e-/hole pairs for  
12.7 keV photon)


# PHOTON COUNTING DETECTORS



Integrated current at  
anode  $\sim$  charge




# PHOTON COUNTING DETECTORS



Total charge  $>$  threshold  
(typically 0.5 photon)

# PHOTON COUNTING BACKGROUND


- Thresholding eliminates read-out noise
- Background caused by atoms in beam which are not crystal
- Well mounted samples will have lower background
- Poorly mounted, or in situ, may have high background



# DETECTOR PARAMETERS

- Pilatus: 172 $\mu\text{m}$  square pixels, 320 $\mu\text{m}$  to 1mm thick
- Eiger: 75 $\mu\text{m}$  square pixels, typically 450 $\mu\text{m}$  thick
- Material: Si common, good at lower energies, transparent at high energy, CdTe becoming more common, excellent at high energy ( $> 20 \text{ keV}$ )
- Pixels are very "small and deep" so spots can smear particularly at higher energies as photons pass through the sensor material

# PHOTON COUNTING DETECTORS - "PILEUP"



Many photons arrive  
within detection interval




# PHOTON COUNTING DETECTORS - "PILEUP"





Much ionisation

# PHOTON COUNTING DETECTORS - "PILEUP"





Lots of charge - but still  
only passes above  
threshold so 1 count

# PHOTON COUNTING DETECTORS - "PILEUP"



Incoming vs. uncorrected count rate  
(photons / s / pixel)

# PHOTON COUNTING DETECTORS - "PILEUP"





Fraction recorded (uncorrected) vs. pixel count in 4ms



# KEY POINTS

- This *is* corrected for but that correction assumes constant count rate in pixel
- If you have finely sliced data, this assumption is probably sound
- If you have recorded data with rotation  $\gg$  mosaic spread the assumption is probably rather poor (more on this later)



# PHOTON COUNTING DETECTORS - EFFICIENCY



i03 Eiger at 10keV

# PHOTON COUNTING DETECTORS - EFFICIENCY

Path length through  
sensor longer  
around edges of  
detector than middle



i03 Eiger at 10keV

# IDEA #1

- Be aware of saturating the detector
- It is quite hard, but a good crystal on a bright beamline is able to do this
- It is hard to see if this has happened
- For a photon counting detector, too weak is rarely a problem
- Be aware of detector efficiency at higher energies
- Talking to beamline staff a good way to go

# DATA INTENSITY / MULTIPLICITY

- The “quality” of data depends on many factors - including the individual spot intensity and number of observations
- Measuring twice as many copies of data half as strong “averages out” to give essentially the same data - with a photon counting detector
- Twice as many, half as strong, measurements also gives opportunity to truncate the data if damage is present



18-11-2021 11:13:02 - gw/20211118/TestInsulin/ins\_1/ins\_1\_1\_master.h5

Sample: ins\_1

Flux: 5.23e+11

Ω Start: 0.0°

Ω Osc: 0.10°

Ω Overlap: 0°

No. Images: 3600

Resolution: 1.51Å

Wavelength: 1.2398Å

Exposure: 0.004s

Transmission: 4.98%

Beamsize: 80x20μm

Type: SAD

Comment: (-303,175,-202)

Aperture: Large

225

200

175

150

125

100

75

0

100

200

300

1.8

2.0

2.2

2.4

2.6

2.8

3.0

Spots

Bragg

Res

Auto Processing

xia2 dials: ✓ xia2 3dii: ✓ fast\_dp: ✓ autoPROC: ✓ autoPROC+STARANISO: ✓

| Type               | Resolution   | Spacegroup | Mn<l/sig(l)> | Rmeas Inner | Rmeas Outer | Completeness | Cell                                | Status                |
|--------------------|--------------|------------|--------------|-------------|-------------|--------------|-------------------------------------|-----------------------|
| xia2 dials         | 54.88 - 1.43 | I 2 3      | 18.1         | 0.037       | 4.378       | 100.0        | 77.62 77.62 77.62 90.00 90.00 90.00 | processing successful |
| xia2 3dii          | 38.85 - 1.51 | I 2 3      | 25.7         | 0.031       | 4.697       | 100.0        | 77.69 77.69 77.69 90.00 90.00 90.00 | processing successful |
| fast_dp            | 27.47 - 1.81 | I 2 3      | 40.0         | 0.029       | 0.827       | 99.5         | 77.69 77.69 77.69 90.00 90.00 90.00 | processing successful |
| autoPROC           | 54.95 - 1.51 | I 2 3      | 25.4         | 0.030       | 5.385       | 100.0        | 77.71 77.71 77.71 90.00 90.00 90.00 | processing successful |
| autoPROC+STARANISO | 54.95 - 1.53 | I 2 3      | 27.7         | 0.029       | 3.182       | 94.7         | 77.71 77.71 77.71 90.00 90.00 90.00 | processing successful |

xia2 dials

xia2 3dii

fast\_dp ⚠

autoPROC

autoPROC+STARANISO

9 check(s) passed

Beam Centre

X

Y

Start

157.61

167.23

Refined

157.56

167.24

Δ

0.05

-0.01

Space Group

A

B

C

α

β

γ

I 2 3

77.62

77.62

77.62

90.00

90.00

90.00

Shell

Observations

Unique

Resolution

Rmeas

l/sig(l)

CC Half

Completeness

Multiplicity

Anom Completeness

Anom Multiplicity

CC Anom

outerShell

17869

722

1.43 - 1.46

4.378

0.3

0.3

100.0

24.7

100.0

12.7

0.0

innerShell

30056

780

3.89 - 54.92

0.037

101.8

1.0

100.0

38.5

100.0

21.4

0.3

overall

552684

14506

1.43 - 54.88

0.103

18.1

1.0

100.0

38.1

100.0

19.8

0.1

Downstream Processing

fast\_ep: ✓ MrBUMP: 2x ✓

x4  
/4

18-11-2021 11:17:40 - gw/20211118/TestInsulin/ins\_1/ins\_1\_3\_master.h5

Sample: ins\_1

Flux: 1.27e+11

Ω Start: -720.0°

Ω Osc: 0.10°

Ω Overlap: 0°

No. Images: 14400

Resolution: 1.51Å

Wavelength: 1.2399Å

Exposure: 0.004s

Transmission: 1.25%

Beamsize: 80x20μm

Type: SAD

Comment: (-303,175,-202)

Aperture: Large

120

100

80

60

40

20

-500

0

500

2.0

2.5

3.0

3.5

4.0

Spots

Bragg

Res

Auto Processing

xia2 dials: ✓ fast\_dp: ✓ xia2 3dii: ✓ autoPROC: ✓ autoPROC+STARANISO: ✓ xia2.multiplex: ✓

| Type               | Resolution   | Spacegroup | Mn<l/sig(l)> | Rmeas Inner | Rmeas Outer | Completeness | Cell                                | Status                |
|--------------------|--------------|------------|--------------|-------------|-------------|--------------|-------------------------------------|-----------------------|
| xia2 dials         | 54.95 - 1.45 | I 2 3      | 21.7         | 0.051       | 10.068      | 100.0        | 77.71 77.71 77.71 90.00 90.00 90.00 | processing successful |
| fast_dp            | 27.49 - 1.96 | I 2 3      | 64.1         | 0.027       | 0.923       | 99.8         | 77.76 77.76 77.76 90.00 90.00 90.00 | processing successful |
| xia2 3dii          | 38.89 - 1.56 | I 2 3      | 34.5         | 0.038       | 4.269       | 100.0        | 77.78 77.78 77.78 90.00 90.00 90.00 | processing successful |
| autoPROC           | 55.01 - 1.55 | I 2 3      | 33.4         | 0.039       | 10.010      | 100.0        | 77.80 77.80 77.80 90.00 90.00 90.00 | processing successful |
| autoPROC+STARANISO | 55.01 - 1.56 | I 2 3      | 36.1         | 0.038       | 7.600       | 94.6         | 77.80 77.80 77.80 90.00 90.00 90.00 | processing successful |
| 2x xia2.multiplex  | 54.92 - 1.41 | I 2 3      | 27.2         | 0.049       | 5.028       | 100.0        | 77.67 77.67 77.67 90.00 90.00 90.00 | processing successful |

xia2 dials

fast\_dp ⚠

xia2 3dii

autoPROC

autoPROC+STARANISO

2x xia2.multiplex

9 check(s) passed

Beam Centre

X

Y

Start

157.61

167.23

Refined

157.56

167.25

Δ

0.05

-0.02

Space Group

A

B

C

α

β

γ

I 2 3

77.71

77.71

77.71

90.00

90.00

90.00

Shell

Observations

Unique

Resolution

Rmeas

l/sig(l)

CC Half

Completeness

Multiplicity

Anom Completeness

Anom Multiplicity

CC Anom

outerShell

85828

708

1.45 - 1.48

10.068

0.3

0.3

100.0

121.2

100.0

62.2

0.0

innerShell

115895

748

3.94 - 54.99

0.051

142.5

1.0

100.0

154.9

100.0

86.2

0.7

overall

2169166

13982

1.45 - 54.95

0.164

21.7

1.0

100.0

155.1

100.0

80.8

0.4



Downstream Processing

fast\_ep: ✓ MrBUMP: 3x ✓

# TABLE 1 NOT SO HELPFUL

| Overall            |              |                |                 |
|--------------------|--------------|----------------|-----------------|
|                    | Overall      | Low resolution | High resolution |
| Resolution (Å)     | 54.88 - 1.43 | 54.92 - 3.89   | 1.46 - 1.43     |
| Observations       | 552684       | 30056          | 17869           |
| Unique reflections | 14506        | 780            | 722             |
| Multiplicity       | 38.1         | 38.5           | 24.7            |
| Completeness       | 100.00%      | 100.00%        | 100.00%         |
| Mean I/σ(I)        | 18.1         | 101.8          | 0.3             |
| R <sub>merge</sub> | 0.102        | 0.037          | 4.289           |
| R <sub>meas</sub>  | 0.103        | 0.037          | 4.378           |
| R <sub>pim</sub>   | 0.017        | 0.006          | 0.874           |
| CC <sub>½</sub>    | 1.000        | 0.999          | 0.288           |

| Overall            |              |                |                 |
|--------------------|--------------|----------------|-----------------|
|                    | Overall      | Low resolution | High resolution |
| Resolution (Å)     | 54.95 - 1.45 | 54.99 - 3.94   | 1.48 - 1.45     |
| Observations       | 2169166      | 115895         | 85828           |
| Unique reflections | 13982        | 748            | 708             |
| Multiplicity       | 155.1        | 154.9          | 121.2           |
| Completeness       | 100.00%      | 100.00%        | 100.00%         |
| Mean I/σ(I)        | 21.7         | 142.5          | 0.3             |
| R <sub>merge</sub> | 0.164        | 0.051          | 10.026          |
| R <sub>meas</sub>  | 0.164        | 0.051          | 10.068          |
| R <sub>pim</sub>   | 0.013        | 0.004          | 0.911           |
| CC <sub>½</sub>    | 1.000        | 1.000          | 0.324           |



## IDEA #2

- Measure more data with less transmission (at least until you know your samples well)
- Easy way to achieve this is stepped transmission experiment - start off weak, then 4x dose, 4x dose etc. (4 x should increase  $1/\sigma(I)$  by  $\sim 2$ )
- Once you know your sample lifetimes, collect more data that are weak to give room for truncation and reduce impact of count rate on data

# COLLECTION ADVICE

*So I am collecting gently, collecting complete and high multiplicity data, what now?*

- Background: ideally minimise this - mount sample carefully
- Also collect data carefully - for example sample alignment



# SAMPLE MOUNTING


- Use suitable mount - similar size and ideally shape as sample
- Ideally only one sample
- Less liquid means better cooling (less ice) and lower background





# MINIMISE BACKGROUND

- The background is under the spot
- If we keep the rotation narrow, we have less background “before and after” the spot
- Strong spots taken from spot finding - most are rather narrow even with  $0.1^\circ$  images (median / mean  $\sim 2$  images)
- $1^\circ$  vs.  $0.1^\circ$  would have 10x background




Histogram of strong spot widths for insulin crystal with  $0.1^\circ$  / frame - most spots 2 pixels wide



# LOOKING AT IMAGES

- Some people collect 1° images because “always have done”
- Some because they understand the data that way
- Better to use tools properly -> learn to use image viewers





# KNOW YOUR EXPERIMENT

- Collecting data is there to meet experiment goals
- If you know what you are trying to do, you can collect “sub optimal” data which still meets your needs but could be more efficient e.g. wider oscillations
- Keep in mind the count rate correction
- All other things being equal however want rotation width  $<$  half mosaic spread



# IDEA #3

- Know your experiment goals and collect data to achieve them
- If you don't know anything specific, try to collect 360°+ of low dose finely sliced data
- Learn to look at the data properly

# RADIATION DAMAGE

- Typically single biggest problem
- Cannot really be corrected for, or undone - perform experiment carefully
- Usually at workshop have “tricky data sets with added radiation damage”
- The data collection is the last experimental step -> everything later can be repeated etc.
- Where possible avoid radiation damage rather than hoping to correct for it in software

# MULTIPLE SAMPLES

- Assumption is your structure / result is indicative of the "truth"
- Combining data from many samples should get you closer to the truth
- Can also help to reduce impact of radiation damage / other experimental effects





# IDEA #4

- Radiation damage big problem - may need to combine data from multiple samples to achieve goals
- Carefully collecting more than you need allows options later
- Avoiding radiation damage much easier than correcting for it

# CONCLUSIONS

- Ask the beamline staff for advice
- Do not destroy your samples
- Detectors can be saturated - lots of weak data better than a little very strong data
- Know your experimental goals - learn how your samples behave then apply this knowledge
- Where possible, avoid radiation damage and don't be afraid of combining data from many samples

# ACKNOWLEDGEMENTS

- Folks from DECTRIS and elsewhere on writing up how the detectors actually work
- Time on Diamond beamlines to collect example data, characterise detectors
- Samples! Thank you to team at Diamond for the steady supply of good quality test crystals



# CONCLUSIONS

- Ask the beamline staff for advice
- Do not destroy your samples
- Detectors can be saturated - lots of weak data better than a little very strong data
- Know your experimental goals - learn how your samples behave then apply this knowledge
- Where possible, avoid radiation damage and don't be afraid of combining data from many samples