



# VMXm: Micro/nanofocus protein crystallography beamline at Diamond


Anna Warren

[anna.warren@diamond.ac.uk](mailto:anna.warren@diamond.ac.uk)



# Introduction

- Target proteins are getting more complex, leading to smaller and more disordered crystals
- Can design ideal experiment to give optimal data quality



# Current Limits

- Dose tolerance of samples cannot be changed – Henderson/Garman limit fixed
- Reduce dose on sample to measure given data quality:
  - Reduce experimental background
  - Cleaner sample mounting
  - Improve analysis for weak and multicrystal data
  - Record rotation data to improve data quality
  - Visualization of micron and sub-micron crystals
- Take advantage of photoelectron escape

# VMXm Aims

- Improve signal to noise by reducing background:
  - Sample environment under vacuum
  - Crystals mounted with minimal liquid
  - Reduce beamsize to match that of the crystal
- Standard rotation data collection on samples down to 500 nm
  - Alignment without the need for X-ray raster scanning
- Optimise sample alignment, sample cooling and data analysis for micron and sub-micron crystals
- Data collections using minimal amounts of sample

# VMXm Specifications





- 6 – 28 keV energy range
- 0.3 – 10  $\mu\text{m}$  (v) & 0.5 – 5  $\mu\text{m}$  (h)
- Flux @12.4 keV:
  - $\sim 10^{12}$  ph/s in  $0.3 \times 3 \mu\text{m}$  (v x h)
  - $> 10^{11}$  ph/s  $0.3 \times 0.5 \mu\text{m}$  (v x h)
- Interchangeable detectors:
  - Pilatus3 6M Si
  - Eiger2 X 9M CdTe



# Sample Preparation Lab



# Sample Characterisation





# Sample Mounting








# Sample Environment

- *In vacuo* sample environment - reduce X-ray background to a minimum
- *In vacuo* sample cooling to preserve samples
- Standard on-axis optical microscope + SEM for sample visualization
- High stability goniometry - permit rotation data to be measured from micro and nanocrystals

# Beamline







Cryo  
gripper

Fluorescence  
detector

Sample  
holder

OAV port

Beamstop XY

Sample XZ  
stages

Goniometer


SED

SEM


Isolated  
cold-stage  
and gripper



# Standard MX




8	5	11	11	15	9	9	4	7	4	13	6	6
10	9	10	13	12	11	7	7	11	8	5	6	8
6	8	10	8	4	5	9	12	9	13	6	6	13
8	7	13	13	10	7	6	9	10	2	9	12	6
10	10	8	8	2	10	11	10	13	6	15	10	7
4	6	8	3	11	23	46	38	22	9	11	7	9
9	11	7	8	12	43	112	83	20	11	9	7	3
8	10	6	7	11	18	24	23	13	11	12	6	6
9	6	7	3	9	13	9	4	7	10	8	4	8
12	7	9	14	7	13	15	8	8	7	4	12	8
11	6	5	9	11	13	12	9	9	10	4	11	11
10	9	10	4	7	7	5	7	11	6	6	4	10
8	8	9	10	7	7	7	14	11	8	7	7	4




90 counts

# VMXm



0	0	1	0	0	0	1	0	0	0	1	0	1
0	0	0	0	0	0	0	0	0	0	0	0	0
0	0	1	0	0	1	1	0	0	0	0	0	1
0	1	0	0	0	1	1	0	1	0	0	0	0
0	0	1	0	1	3	2	1	0	0	0	0	0
0	0	3	0	1	2	3	1	0	0	0	0	0
0	0	1	1	0	3	2	3	1	0	0	0	0
0	2	1	0	0	3	3	2	0	1	0	0	1
0	0	0	0	1	1	2	0	0	0	0	0	0
0	0	0	0	0	1	1	1	1	0	0	0	0
1	0	1	0	0	0	0	2	0	0	0	0	0
0	0	1	0	1	0	0	0	0	0	1	0	0
0	0	0	0	0	0	0	1	0	0	0	0	0




3 counts


# CPV Ld14 – Cytoplasmic polyhedrosis virus

- Spacegroup I23
- Unit cell  $a=b=c=103 \text{ \AA}$
- Approximate crystal size 2 – 4  $\mu\text{m}$
- Data collected at 21.3keV – exploit photoelectron escape
- Eiger2 X CdTe 9M

	Ji et al., I24	VMXm
Detector	Pilatus 6M Si	Eiger2 X CdTe 9M
Number of crystals	20	1
Energy keV	12.458 keV	21.3 keV
Resolution	72.5–1.91 (2.02–1.91)	72.9–1.30 (1.32–1.30)
Unique Reflection	12952(1045)	40616(1350)
Completeness (%)	92.3(52.1)	90.7(62.2)
$R_{\text{merge}}$	0.199(0.327)	0.284(0.881)
$I/\sigma I$	8.4(1.8)	3.8(0.5)
CC1/2	??	0.933(0.407)
Beamsize	6 x 6 $\mu\text{m}$ slit	3.6 x 3.6 $\mu\text{m}$




Stat	Value
Resolution	72.94–1.38
N. Reflections all/free	37547/1920
R/Rfree	0.136/0.193
RMS dev	
Bonds	0.0141
Angles	1.806






# CPV Us17 – Cytoplasmic polyhedrosis virus

- Spacegroup I23
- Unit cell  $a=b=c=105 \text{ \AA}$
- Approximate crystal size 1-2  $\mu\text{m}$
- Data collected at 21.3keV
- Eiger2 X CdTe 9M
- Merged 17 datasets



	VMXm	XFEL	I24
Number of crystals	13	5787	768
Resolution	74.28–1.46 (1.49–1.46)	28.3–1.75 (1.79–1.75)	74.16–2.20 (2.26–2.20)
Unique Reflection	32555(1664)	19096(1032)	9376(931)
Completeness (%)	97.2(100)	100(100)	99.9(100)
$R_{\text{merge}}$	9.265(- 14.862)	-	0.665(0.000)
$I/\sigma I$	6.7(0.8)	-	6.4(1.4)
CC1/2	0.469(0.268)	-	-









# Mpro ligand co-crystals

- Ligand absent from large crystal form
- Co-crystals with ligand of needle form grown
- Needle crystals (3-4  $\mu\text{m}$  wide)
- 40° wedges of diffraction collected from 12 crystals
  - 3 x 3  $\mu\text{m}$  beam at 21.3 KeV
- dials.multiplex used to combine the 12 datasets

Resolution (Å)	22.53 – 1.81
Wavelength (Å)	0.5821
Observations	371253 (6992)
Unique Observations	23222 (1081)
Multiplicity	16.0 (6.5)
Completeness	81.34% (76.83%)
Mean I/ $\sigma$ I	6.2 (0.6)
Rmeas	0.511(4.039)
R <sub>p</sub> im	0.110 (1.386)
CC1/2	0.98 (0.031)





Stat	Value
Resolution	22.53 – 1.81
N. Reflections all/free	21520/1140
R/R <sub>free</sub>	0.19/0.24
RMS dev	
Bonds	0.01
Angles	1.459




# ToNV


- Polyhedra protein from a nudivirus
- Self-assembles into a dense lattice around new viral particles
- 3 x Met engineered into WT clone to allow SeMet protein
- P3<sub>1</sub>21 a=53.5 Å, c=105.2 Å,  $\gamma=120^\circ$
- Solvent content 21%
- Crystals ~ 5 – 7  $\mu\text{m}$  (SeMet); 3 – 5  $\mu\text{m}$  (WT)



## Se-Met to 1.86Å



## WT to 1.7Å



	Se-Met TonV on VMXm
Detector	Eiger2 X CdTe 9M
Number of crystals	67
Energy keV	12.67 keV
Resolution	105.0-1.86 (1.90-1.86)
Unique Reflection	13841(244)
Completeness (%)	91.3(32.0)
R <sub>merge</sub>	0.288(1.342)
I/σI	17.7(0.7)
CC1/2	0.971(0.133)
Beamsize	3.6 x 3.6 $\mu\text{m}$

# Final Remarks

- VMXm can record high quality diffraction data where both size and number of crystals are limited
- Exploiting low background, sample mounting, high energy and photoelectron escape to optimize experiment
- Successfully demonstrated the VMXm principal of operation
  - Rotation data recorded from crystals a few microns in size
  - Crystals mounted on cryoEM grids in vacuum at cryo-temperatures



# Acknowledgements

Graham Duller  
Andy Stallwood  
Richard Parr  
Tony Gardner  
Dan Magrath  
Marcell Nagy

Leo Cahill  
Richard Littlewood  
Kevin Wilkinson  
Alan Nash  
Steve Buxcey  
Paul Benson's IFM Team

Andy Foster  
James O'Hea  
Brian Nutter  
Ben Kemp

Victoria Lawson  
Neil Smith  
Paul Coe

Mark Lunnon  
Dave Butler  
Adam Taylor  
Adam Prescott  
Leon Adams  
Thomas Hartrampf

Simon Lay  
Hugo Shiers  
Guenther Rehm  
Chris Bloomer  
Sandira Gayadeen

David Laundry  
Lucia Alianelli  
John Sutter  
Simon Alcock  
Kawal Sahwney

Sonia Moon



Gwyndaf  
Evans



Jose  
Trincao



Adam  
Crawshaw

Contact a member of the  
VMXm team if interested to  
know more or use the  
beamline

