


# ProSMART

Injection of prior knowledge to aid new structure determination


- **External Restraints from homologous structures**
  - Protein or nucleic acid chains
- **Hydrogen bond restraints**
  - Protein backbone
- **Generic self-restraints**
  - Everything – protein, nucleic acid, ligand, water
- **Structure analysis**
  - Alignment & comparison - helps analyse differences between models

***Independent of global conformation***

# ProSMART External Restraints



Prior information:




Stabilises structural features


3g4w - 3.7 Å

# External Restraint Generation

structure to be refined



known similar structure (prior)




(abstract representation of an atomic model; circles = atoms)

# External Restraint Generation

structure to be refined

known similar structure (prior)




(abstract representation of an atomic model; circles = atoms)



# External Restraint Generation


structure to be refined



$$d \sim N(r, \sigma^2)$$

(abstract representation of an atomic model; circles = atoms)


# ProSMART Restraints





# Robust Estimation

Cyan: original




# Robust Estimation





# Robust Estimation



Cyan: original  
Green: homolog  
Yellow: refined

Red: long  
Grey: similar  
Blue: short



# Robust Estimation




Cyan: original  
Green: homolog  
Yellow: refined

Red: long  
Grey: similar  
Blue: short

# Motivational Example

Ovotransferrin



1ryx - 3.5Å

Low-resolution refinement:

Weak signal  
Noisy data




Unstable refinement

**Result:**  
Poor quality model


# Motivational Example

Ovotransferrin



1ryx - 3.5Å

High-resolution homologue




2d3i - 2.15Å




# Motivational Example

Ovotransferrin



Models don't superpose well

# Example: Ovotransferrin




1ryx (3.5Å)

# Example: Ovotransferrin

Restraints:  
Backbone  
Side chains



1ryx (3.5Å)  
restrained to  
2d3i (2.15Å)



Red: long  
Grey: similar  
Blue: short


# Example: Ovotransferrin

Restraints:  
Backbone  
Side chains



After re-refinement


1ryx (3.5Å)  
restrained to  
2d3i (2.15Å)



Red: long  
Grey: similar  
Blue: short


# External Restraints

Ovotransferrin



# External Restraints


Ovotransferrin





# External Restraints

Ovotransferrin




# External Restraints

Ovotransferrin

Original Structure

R/R<sub>free</sub> : 0.286/0.330






# External Restraints

Ovotransferrin


Original Structure

$R/R_{\text{free}}$  : 0.286/0.330




Re-refined with External Restraints

$R/R_{\text{free}}$  : 0.263/0.307



# External Restraints



1.3 $\sigma$

Original Structure

R/R<sub>free</sub> : 0.286/0.330




External restraints

(40 cycles)

R/R<sub>free</sub> : 0.263/0.307



# External Restraints



1.3 $\sigma$

Original Structure

R/R<sub>free</sub> : 0.286/0.330



External restraints

(40 cycles)

R/R<sub>free</sub> : 0.263/0.307



Modify

Real Space Refine




Jelly body

(40 cycles)

R/R<sub>free</sub> : 0.253/0.304



# External Restraints



1.3 $\sigma$

Original Structure

R/R<sub>free</sub> : 0.286/0.330




External restraints

(40 cycles)

R/R<sub>free</sub> : 0.263/0.307



# External Restraints



1.3 $\sigma$

Original Structure

R/R<sub>free</sub> : 0.286/0.330



External restraints

(40 cycles)

R/R<sub>free</sub> : 0.263/0.307



Build TYR92

Modify LYS209





Jelly body

(40 cycles)



R/R<sub>free</sub> : 0.252/0.307

# LibG Nucleic Acid Restraints

Base-pair restraints




Parallel plane restraints





# LibG Nucleic Acid Restraints



Red: long  
Grey: similar  
Blue: short

Visualise & edit in Coot

# External Restraints

When refining at low resolution, check:

- Refinement statistics – *Not always conclusive*
- Geometry – *Not always conclusive*
- Electron density – *Not always reliable*

**Conclusion:** At low resolution, everything has to add up!  
Take care; reflect

Quality of prior information is important – consider manual re-refinement  
– PDB-REDO is useful





# Automated pipeline - LORESTR

- Efficiency of ProSMART-generated restraints greatly depends on the homologues used
- If several homologues are available, substantial manual effort is required to find their optimal combination
- Other refinement parameters (scaling, solvent, etc) also affect efficiency of the process


Solution:

**LOW-Resolution Structure Refinement**


# Automated pipeline - LORESTR

**Low-resolution Structure Refinement Pipeline**


# Automated pipeline - LORESTR




# Automated pipeline - LORESTR




# Automated pipeline - LORESTR



# Automated pipeline - LORESTR




# Automated pipeline - LORESTR







# Automated pipeline - LORESTR




# Automated pipeline - LORESTR




# Automated pipeline - LORESTR




# Automated pipeline - LORESTR



# Automated pipeline - LORESTR



# Automated pipeline - LORESTR



# Summary

## Tools to help with model building and refinement:

**REFMAC5:** Refinement, jelly body restraints, map sharpening/blurring

**ProSMART:** External restraints, comparative analysis

**LibG:** Nucleic acid restraints

**LORESTR:** Automated low-resolution pipeline

**AceDRG:** Ligand dictionary and conformer generation

**Coot:** Visualisation & manipulation of restraints, map blurring  
...also morphing, jiggle-fit, backrub rotamers...

*Many tools are applicable to cryo-EM as well as MX*

# What and When

## **Early stages (e.g. straight after MR)**

- Rigid body refinement
- Jelly body – sometimes up to 200 cycles

## **Medium stages – during model building**

- Auto local NCS – wherever possible
- External restraints (40 cycles) – homologue available
- Otherwise, jelly body... but not together
- H-bond and DNA/RNA restraints – no homologue available
- Secondary structure conformation restraints – model building tool
- Add hydrogens (?)

## **Medium-final stages**

- TLS – at medium resolutions
- Anisotropic B-factors – only at high resolution
- Twin refinement – only if you are sure

## **Final stages of refinement**

Jelly body – around 20 cycles



# Relevant Publications

## Low-resolution refinement with REFMAC5, ProSMART, LibG & LORESTR:

- Nicholls *et al.* (2017) ) Low Resolution Refinement of Atomic Models Against Crystallographic Data. *Protein Crystallography*, 565-93.
- Nicholls *et al.* (2013) Recent Advances in Low Resolution Refinement Tools in REFMAC5. *Adv. Methods for Bio. Xtallography*, 231-58.
- Nicholls *et al.* (2012) Low Resolution Refinement Tools in REFMAC5. *Acta Cryst.* D68, 404-17.

## Tools for cryo-EM model fitting & refinement:

- Nicholls *et al.* (2018) Current approaches for the fitting and refinement of atomic models into cryo-EM maps using CCP-EM. *Acta Cryst.* D74, 492-505.
- Murshudov (2016) Refinement of atomic structures against cryo-EM maps. *Methods in Enzymology*, 277-305.
- Brown *et al.* (2015) Tools for macromolecular model building and refinement into electron cryo-microscopy reconstructions. *Acta Cryst.* D71, 136-53.

## Tools for ligand fitting & validation:

- Nicholls (2017) Ligand fitting with CCP4. *Acta Cryst.* D73, 158-170.
- Emsley (2017) Tools for ligand validation in Coot. *Acta Cryst.* D73, 203-10.
- Debreczeni & Emsley (2012) Handling ligands with Coot. *Acta Cryst.* D68, 425-30.

## Cooperative utilisation of information from Xtal and NMR:

- Kovalevskiy *et al.* (2018) Overview of refinement procedures within REFMAC5: Utilising Data from Different Sources. *Acta Cryst.* D74, 215-27.
- Carlon *et al.* (2016) How to tackle protein structural data from solution and solid state: An integrated approach. *Progress in nuclear magnetic resonance spectroscopy*. 92, 54-70.


## Effect of Twinning on R-factors:

- Murshudov GN (2011) Some properties of Crystallographic Reliability Index – Rfactor: Effect of Twinning. *Appl. & Comp. Math.*, 10, 250-61.

# Acknowledgements

Contact: [lucrezia@mrc-lmb.cam.ac.uk](mailto:lucrezia@mrc-lmb.cam.ac.uk)  
[www2.mrc-lmb.cam.ac.uk/groups/murshudov/](http://www2.mrc-lmb.cam.ac.uk/groups/murshudov/)

## *MRC-LMB Computational Structural Biology Group*



Rob Nicholls  
Fei Long  
Oleg Kovalevskiy  
**Garib Murshudov**  
Rangana Warshamanage  
Paul Emsley  
Keitaro Yamahita

### ***CCP4 Core***

Eugene Krissinel  
Andrey Lebedev  
Charles Ballard  
Ronan Keegan  
Ville Uski

### ***CCP4i2***

Martin Noble  
Stuart McNicholas  
Jon Agirre  
Liz Potterton

### ***CCP-EM***

Martyn Wynn  
Tom Burnley  
Colin Palmer

### ***Global Phasing***

Marcin Wojdyr

### ***Collaborators***

Marcus Fischer  
Robbie Joosten  
Andrea Thorn  
**Roberto Steiner**  
Alan Brown  
Jude Short  
Ana Casañal  
Rafiga Masmaliyeva  
Azzurra Carlon

### ***Computing***

Jake Grimmett  
Toby Darling

All colleagues from  
MRC-LMB, CCP4, CCP-EM  
Users for feedback!