

DLS-CCP4 Data Collection and Structure Solution Workshop
6 November 2023

MX at Diamond

Marco Mazzorana



Diamond and the Harwell Campus



What is Diamond

The UK national synchrotron
(86% STFC, 14% Wellcome Trust)

20 years from foundation

15 years of operations



What is Diamond

Phase 1
7 beamlines
completed 2007

Phase 2
+15 beamlines
completed 2012

Phase 3
+11 beamlines to
be available by
2020



Many different instruments...

33 beamlines

8 Science Groups

2 national cryo-EM facilities


12k+ users

100+ companies


Structural biology labs

MPL, crystallization facility XChem


UK-XFEL Hub



...to investigate the structure of biological systems



The power of combining techniques



Super resolution
cryo-light
microscopy

Cryo X-ray
tomography

Cryo-electron
tomography


Atomic models
from MX fit to
averaged particles
of two
intermediates

Understand interactions and action across a large range of spatial and time domains (often weak and transient)

Reveal inner workings of the cell with a full dynamic picture of living organisms

Molecular description of cellular organization can transform fundamental science and understanding of disease processes

What crystallography tells us about macromolecules



Challenges and synergies

- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus



- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long Wavelength



- VMXi
- I24
- I03
- Screening
- Data collection
- CLIII

In situ



- Long wavelength
- Multi-axis strategies
- Inverse beam
- Wedged MAD
- Pipelines

Experimental Phasing



- *In situ*
- Containment Level II & III
- Schedule 5

Biological Containment



- Humidity control
- Spectroscopy
- Tomography
- Dynamics
- Fluorescence

Complementary information



- Diffraction image processing for SR Electron and XFEL data
- Collaborative effort
- Led and managed by DLS

Dials Project




- Experiment database
- Synchweb
- Data archive
- Data access
- Data reprocess

Data Management




Dealing with the smallest crystals: microfocus


3 microfocus beamlines



I04



I24



VMXm

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*


Microfocus




Microfocus for time-resolved studies

Sample delivery for serial crystallography (SSX and SFX)


Liquid jets/ LCP extruders



Acoustic droplet ejection



Fixed target




- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




Microfocus for time-resolved studies

Fixed-target at I24




Bespoke mount in I24
Fast stages



Full chip rastered in 6-20 mins
10-50 ms/frame (+ pause at each frame)
Typically ~25 chips/12 hours

Multi-crystal structures (from stills)



1.5 Å structure
from 9000 crystals

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus



The UK-XFEL Hub



POHANG ACCELERATOR LABORATORY



- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




tools and strategies to enable dynamic structural biology facilitating:

- (time-resolved) serial structural biology experiments via sample preparation, delivery, data collection, and processing
- the transfer of methods between XFEL, synchrotron, and/or cryo-EM sources
- access to, and data collection from complementary facilities (DLS, CLF, eBIC)



Extending the time range




- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




The Diamond-II project

4th generation synchrotrons and XFELs have hotter beams




Smaller than visible

VMXm: the first nano-focus beamline



< 500 nm beamsize
Integrated SEM
In vacuum
Tunable 7-25 keV (2.0-0.57Å)



- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




Smaller than visible


VMXm: the first nano-focus beamline

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*


Microfocus




EM-grids as sample holders



OAV




SEM



extremely low background

Smaller than visible

VMXm: the first nano-focus beamline




27 crystals

1.5-2.0 deg/data per crystal

Rigid body+restrained refinement

Dials, CCP4, Refmac, Coot



Summary of merging statistics

	Overall	Low resolution	High resolution
Resolution (Å)	36.57 - 1.50	36.58 - 4.07	1.53 - 1.50
Observations	116543	7779	1406
Unique reflections	25790	1405	833
Multiplicity	4.5	5.5	1.7
Completeness	87.23%	90.30%	57.21%
Mean I/ σ (I)	5.1	13.4	0.7
Rmerge	0.105	0.128	0.775
Rmeas	0.217	0.140	1.017
Rpim	0.088	0.053	0.646
CC1/2	0.969	0.970	0.319

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




The need for longer wavelength

Biologically relevant elements

	e (keV)	$\lambda(\text{\AA})$
P	2.146	5.779
S	2.472	5.016
Cl	2.822	4.393
K	3.607	3.437
Ca	4.038	3.070

Xray absorption edges outside the capabilities of standard beamlines




- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long Wavelength



Challenges and solutions



Long wavelength increases the angle of diffraction
Need a curved detector to measure high resolution spots

- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems


Long
Wavelength



Challenges and solutions

Absorption by air of
long wavelength X-rays

In vacuum beamline



Bespoke sample holders




- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems


Long
Wavelength




Results: metal identification



Peptidyl transferase center



Decoding center




- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems


Long
Wavelength




Extending the element 'palette'




ribosome 70S




Na/K transporter




FAD-binding



SLC7 transporter



PETase



- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long
Wavelength



High throughput



BART sample changer:
37 unipucks
597 samples
Exchange rate approx. 20 sec/sample

EIGER2 XE 16M

Pixel size: $75 \times 75 \mu\text{m}^2$

18,093,576 pixel

560 frames/sec (3600 frames in <7 sec)



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



Logistics



Weland/Randex automated store

high capacity (up to 272 dewars)
barcode tracking (dewars and pucks)
logistics integrated in Synchweb/ISPyB




- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



Responsive/dynamic scheduling

2019



2020



2022




- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



Unattended data collection

Optimized for speed




- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)


High Throughput



Xray-centering and strategies



Automated processing pipelines



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)



High Throughput



HT fragment screening



Conventional screening

Preselection of best compounds



Fragment-based approach

Screen by crystal



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)


High Throughput




HT fragment screening

XChem Lab at Diamond


Fragment libraries




Beamline I04-1



~700 datasets/day, unattended



Protein crystals in plates



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



1 crystal – 1 fragment

Soak/harvest/collect data
(up to 1000 crystals/week)

Routine since 2016

March 2020 COVID-19 targets

Access to Diamond



Non-proprietary (academic)

Peer-review access based on quality of proposal

- Free access
- Publishable experiments

Access modes:

- **BAG** (2 year programme)
- **Single application** (6 months)
- **Rapid access** (apply anytime for one-off experiments)

MX beamlines, XChem



Proprietary (industry and academia)

Organized via Industrial Liaison Office

- No need to publish
- Allocated slots throughout the schedule
- 10% of allocated beamtime

Services provided:

- Direct and remote access
- Unattended data collection
- Mail-in data collection and analysis service
- XChem fragment screening



The team

I04-1 / XChem

- Frank von Delft
- Jose Brandao-Neto
- Louise Dunnett
- Kutumbarao Nidamarthi
- Daren Fearon
- Charlie Thomlinson

I03

- Dave Hall
- Neil Paterson
- Mark Williams
- Felicity Bertram

I04

- Ralf Flaig
- David Aragão
- Marco Mazzorana
- Kutumbarao Nidamarthi
- Pierpaolo Romano

VMXm

- Gwynndaf Evans
- Jose Trincão
- Anna Warren
- Adam Crawshaw

VMXi

- Mike Hough
- Juan Sanchez-Weatherby
- James Sandy
- Halina Mikolajek

XFEL-Hub

- Allen Orville
- Pierre Aller
- Anastasiia Shilova
- Tiankun Zhou
- Jos Kamps

I24

- Robin Owen
- Danny Axford
- Sofia Jaho
- Do-Heon Gu
- Pierre Aller
- Anastasiia Shilova

I23

- Armin Wagner
- Kamel el Omari
- Vitaliy Mykhaylyk
- Ramona Duman
- Christian Orr

MPL

- Andrew Quigley
- James Birch
- Harish Cheruvuvara

MX User Support

- Marco Mazzorana
- Kutumbarao Nidamarthi
- Felicity Bertram

Data Analysis & Management

- Graeme Winter
- Irakli Sikharulidze
- Neil Smith
- Karl Levik
- Stu Fisher

DIALS

- Nicholas Devenish
- James Beilsten-Edmonds
- Luis Fuentes-Montero
- Graeme Winter
- David Waterman

DIALS
Diffraction Integration for Advanced Light Sources