



# Getting the most from micron sized crystals


Anna Warren

[anna.warren@diamond.ac.uk](mailto:anna.warren@diamond.ac.uk)










# Introduction

- X-ray crystallography remains the gold standard for the determination of macromolecular structure and protein substrate interactions
- Production of suitable crystals remains one of the largest bottlenecks
- Many crystallisation strategies have been developed
- Synchrotron beamlines have also improved to assist biochemists with their structural problems
- However, some complex proteins are difficult to crystallise and don't form large uniform crystals



# Introduction

- Sometimes only small crystals are formed
- Can design ideal experiment to give optimal data quality



# Previous Limits

- 2.2 Å data can be collected from 1  $\mu\text{m}^3$  crystals (~700 well diffracting crystals)
- 3 Å data can be collected from 5  $\mu\text{m}^3$  membrane protein crystals grown in LCP (~35 crystals, grid scanned first for centring)
- From theoretical calculations a complete 2 Å dataset can be collected from a single 1  $\mu\text{m}^3$  lysozyme crystal (Holton and Frankel, 2010)
- Discrepancies between theory and experiment




# Current Limits

- Dose tolerance of samples cannot be changed – Henderson/Garman limit fixed
- Reduce dose on sample to measure given data quality:
  - Reduce experimental background
  - Cleaner sample mounting
  - Improve analysis for weak and multicrystal data
  - Record rotation data to improve data quality
  - Visualization of micron and sub-micron crystals
- Take advantage of photoelectron escape

# VMXm Aims

- Improve signal to noise by reducing background:
  - Sample environment under vacuum
  - Crystals mounted with minimal liquid
  - Reduce beamsize to match that of the crystal
- Standard rotation data collection on samples down to 500 nm
  - Alignment without the need for X-ray raster scanning
- Optimise sample alignment, sample cooling and data analysis for micron and sub-micron crystals
- Data collections using minimal amounts of sample





# Sample Preparation




# Sample Preparation Lab











# Cryo-EM Grids

3 mm





Support film









# Sample Blotting





# Sample Characterisation





# Sample Mounting








# Beamline





# Beamline





Cryo  
gripper

Fluorescence  
detector

Sample  
holder

OAV port

Beamstop XY

Isolated  
cold-stage  
and gripper

Sample XZ  
stages

Goniometer

SED

SEM









# Beamline Results



# Standard in air MX





# VMXm in vacuo




# CPV Us17 – Cytoplasmic polyhedrosis virus

Previously measured on I24 and at CXI instrument @LCLS. Initially solved my MR using LCLS data to 1.75 Å and then extended with better XFEL data analysis to 1.46 Å.

- Cytoplasmic polyhedrosis virus
- Spacegroup I23
- Unit cell  $a=b=c=105$  Å
- Approximate crystal size 1-2  $\mu\text{m}$
- Data collected at 21.3keV
- Eiger2 X CdTe 9M





	VMXm	XFEL	I24
Number of crystals	14	6537	768
Resolution	74.2–1.45 (1.48-1.45)	25.0-1.46 (1.485-1.46)	74.16-2.20 (2.26-2.20)
Unique Reflection	33145(1699)	34369 (-)	9376(931)
Completeness (%)	97.2(100)	99.5(92.8)	99.9(100)
$R_{\text{meas}}$ (%)	89.1(668.1)	19.0(71.4)	0.665(0.000)*
$I/\sigma I$	9.2(1.1)	3.2 (-)	6.4(1.4)
CC1/2	0.986(0.302)	0.999(0.331)	--
$R_{\text{work}}/R_{\text{free}}$ (%)	15.2/16.6	11.1/15.8	14.7/19.9






# ToNV


- Polyhedra protein from a nudivirus
- Self-assembles into a dense lattice around new viral particles
- 3 x Met engineered into WT clone to allow SeMet protein
- $P3_121$   $a=53.5$  Å,  $c=105.2$  Å,  $\gamma=120^\circ$
- Solvent content 21%
- Crystals  $\sim 5 - 7$   $\mu\text{m}$  (SeMet);  $3 - 5$   $\mu\text{m}$  (WT)



## Se-Met to 1.86Å



## WT to 1.7Å






	Se-Met TonV on VMXm
Detector	Eiger2 X CdTe 9M
Number of crystals	67
Energy keV	12.67 keV
Resolution	105.0-1.86 (1.90-1.86)
Unique Reflection	13841(244)
Completeness (%)	91.3(32.0)
$R_{\text{merge}}$	0.288(1.342)
$I/\sigma I$	17.7(0.7)
CC1/2	0.971(0.133)
Beamsize	3.6 x 3.6 $\mu\text{m}$




# Mpro ligand co-crystals

- Ligand absent from large crystal form
- Co-crystals with ligand of needle form grown
- Needle crystals (3-4  $\mu\text{m}$  wide)
- 40° wedges of diffraction collected from 12 crystals
  - 3 x 3  $\mu\text{m}$  beam at 21.3 keV
- dials.multiplex used to combine the 12 datasets

Resolution (Å)	22.53 – 1.81
Wavelength (Å)	0.5821
Observations	371253 (6992)
Unique Observations	23222 (1081)
Multiplicity	16.0 (6.5)
Completeness	81.34% (76.83%)
Mean I/ $\sigma$ I	6.2 (0.6)
R <sub>meas</sub>	0.511(4.039)
R <sub>pim</sub>	0.110 (1.386)
CC1/2	0.98 (0.031)




Stat	Value
Resolution	22.53 – 1.81
N. Reflections all/free	21520/1140
R/R <sub>free</sub>	0.19/0.24
RMS dev	
Bonds	0.01
Angles	1.459




# Membrane protein - FFAR1

- GPCR membrane protein – crystals  $\sim 5\ \mu\text{m}$
- Grown in classic glass plate sandwich in LCP – highly viscous, adds significant background scatter
- Crystallisation solution including 12% MPD used to wash the bolus out from under the glass
- Applied directly to grid mounted in plunge freezer as per soluble sample



	Overall	Low resolution	High resolution
Resolution ( $\text{\AA}$ )	53.79 - 2.27	53.81 - 6.16	2.31 - 2.27
Observations	297971	14118	15053
Unique reflections	30700	1593	1512
Multiplicity	9.7	8.9	10.0
Completeness	99.96%	99.44%	100.00%
Mean $I/\sigma(I)$	6.3	26.3	0.6
$R_{\text{merge}}$	0.285	0.106	4.440
$R_{\text{meas}}$	0.302	0.113	4.687
$R_{\text{pim}}$	0.096	0.036	1.455
$CC_{1/2}$	0.990	0.993	0.277





# Current status

- VMXm can record high quality diffraction data where both size and number of crystals are limited
- Exploiting low background, sample mounting, high energy and photoelectron escape to optimize experiment

## User programme:

- Current commissioning call will continue until further notice
- Sample size envelope is broad – not just micron sized
  - Even though beam size is smaller we can sensibly prepare grids for the beamline with larger crystals
- Send VMXm staff pictures of sample!
  - This will help us advise on suitability and approach
- Drops with low numbers of crystals are OK! Do not need high concentration of crystals!
  - Access is relatively light touch at present requiring a brief scientific justification

# VMXm Team



Gwyndaf  
Evans



Jose  
Trincao



Adam  
Crawshaw

Contact a member of the VMXm team if interested to know more or use the beamline

**For the DLS-CCP4 Workshop data collection day:  
Monday 27<sup>th</sup> November**

[Anna.Warren@diamond.ac.uk](mailto:Anna.Warren@diamond.ac.uk)  
[VMXM@diamond.ac.uk](mailto:VMXM@diamond.ac.uk)

