


XDS usage - introduction and demonstration

Kay Diederichs



Outline

- 1) Introduction to XDS (XSCALE, XDSCC12)
- 2) Demo: Processing data with XDSGUI (and using XDSSTAT, XDSCC12, SPOT2PDB, COOT, POINTLESS)

The *XDS* program suite

- Original author: Wolfgang Kabsch (Max-Planck-Institute Heidelberg)
- Since ~1986
- I joined 2007



The XDS+ programs

- **XDS**: the main program - indexing, integrating, scaling, statistics
- **XSCALE**: scale several XDS intensity data sets together; zero-dose extrapolation; statistics
- **XDSCONV**: convert to MTZ / SHELX /... format (AIMLESS and CTRUNCATE are not needed!)

Programs independent from the XDS distribution:

- **XDS-Viewer**: inspect diagnostic images written by XDS, or (single) data frames (open source). *adxv* or *dials.image_viewer* may be used instead.
- **XDSTAT**: additional statistics
- **XDSGUI**: graphical user interface for XDS, SHELX C/D/E, ARCIMBOLDO (open source)
- **XDSCC12**: (XDS) which frames are bad?
(XSCALE) which data sets to re-index and merge?
- **XSCALE_ISOCLUSTER** multi-data-set: visualize relations and cluster

Sources of information

- XDS main website: <https://xds.mr.mpg.de> ; complete, accurate, up-to-date documentation; download: Linux, Intel- and Silicon-Mac; Windows - use WSL.
- XDWiki: <https://wiki.uni-konstanz.de/xds>
Installation; data sets; documentation; download; links to e.g. Matthew J. Whitley's excellent tutorial given at CSHL 2018
- CCP4 bulletin board
- SBGrid talk (May 2020) at <https://www.youtube.com/watch?v=3WU9NrILECo>
- Making a difference in multi-data-set crystallography: simple and deterministic data-scaling/selection methods. Assmann, G.M., Wang, M., Diederichs, K. (2020) Acta Cryst D76, 636 (serial crystallography, XDSCC12)

Automatic processing with XDS

- beamline software (provides **XDS.INP**)
- scripts: **xia2** (CCP4), **autoPROC** (Globalphasing), **generate_XDS.INP** (XDSwiki), **fast_dp** (Diamond), *xdsme* (Soleil), *autoxds* (SSRL), *autoprocess* (CMCF), ...
- CCP4: *pointless*, *xdsconv* (type CCP4_I+F, or CCP4, or CCP4_I, or CCP4_F)
- SHELX: *shelxc* reads XDS_ASCII.HKL

Principle of XDS processing

- There is one JOB= line in **XDS.INP** which specifies a list of tasks:

JOB= XYCORR INIT COLSPOT IDXREF DEFPIX INTEGRATE CORRECT


- data reduction is divided into tasks in a **modular** way
- information storage/exchange/flow between tasks by data files which may be inspected/analyzed
- each task needs the result from the previous tasks
- fine-tuning of a task does *not* require previous tasks to be repeated
- each task writes its output file **<TASK>.LP**

The tasks are ...

- XYCORR : write positional correction files
(**X-CORRECTIONS.cbf**, **Y-CORRECTIONS.cbf**)
- INIT : find background pixels (defaults usually OK)
- COLSPOT: find reflection positions
- IDXREF : "index" reflections; user may supply/choose spacegroup
- XPLAN [not required] : strategy for data collection
- DEFPIX : mask shadows on detector (use *XDSGUI!*)
- INTEGRATE : evaluates intensities on all frames, writes **INTEGRATE.HKL** and **FRAME.cbf**
- CORRECT : **scales**, rejects outliers, statistics, writes scaled, unmerged **XDS_ASCII.HKL** (and other files)

Information flow

NAME_TEMPLA	OSCILLATION_	ORGX	
TE_OF_DATA_	RANGE	ORGY	
FRAMES	SEPMIN	DETECTOR_DISTANCE	DATA_RANGE
DETECTOR	STRONG_PIXEL	X_RAY_WAVELENGTH	
		SPACE_GROUP_NUMBER	



Example **XDS.INP**

```
JOB= XYCORR INIT COLSPOT IDXREF DEFPIX INTEGRATE CORRECT
ORGX=1546 ORGY=1552      !Detector origin (pixels); e.g. NX/2 NY/2
DETECTOR_DISTANCE=180    ! (mm)
OSCILLATION_RANGE=0.50   !degrees (>0)
X-RAY_WAVELENGTH=0.980243 !Angstroem
NAME_TEMPLATE_OF_DATA_FRAMES=frms/wga2-27_1_???.img      !No blanks allowed!
DATA_RANGE=1 360      !Numbers of first and last data image collected
BACKGROUND_RANGE=1 10 !Numbers of first and last data image for background
SPACE_GROUP_NUMBER= 19 !0 for unknown crystals; cell constants are ignored.
UNIT_CELL_CONSTANTS= 44.4 86.4 104.5 90 90 90 ! not required if spgr=0
REFINE(IDXREF)=BEAM AXIS ORIENTATION CELL DISTANCE
REFINE(INTEGRATE)=DISTANCE BEAM ORIENTATION CELL ! AXIS
ROTATION_AXIS= 1.0 0.0 0.0
INCIDENT_BEAM_DIRECTION=0.0 0.0 1.0
FRACTION_OF_POLARIZATION=0.99      ! SLS X06SA
POLARIZATION_PLANE_NORMAL= 0.0 1.0 0.0
DETECTOR=CCDCHESS      MINIMUM_VALID_PIXEL_VALUE=1      OVERLOAD=65000
DIRECTION_OF_DETECTOR_X-AXIS= 1.0 0.0 0.0
DIRECTION_OF_DETECTOR_Y-AXIS= 0.0 1.0 0.0
VALUE_RANGE_FOR_TRUSTED_DETECTOR_PIXELS= 7000 30000 !Used by DEFPIX
      !for excluding shaded parts of the detector.
INCLUDE_RESOLUTION_RANGE=50.0 1.3 !Angstroem; used by DEFPIX,INTEGRATE,CORRECT
```

Bold keyword/parameter pairs are required. Complete documentation at
xds.mr.mpg.de/html_doc/xds_parameters.html

Example **XSCALE.INP**

```
!===== EXAMPLE 3: specific reindexing of input data sets
!
!      Use of specific reindexing of input data sets for resolving
!      indexing ambiguities in the scaled output data set. This
!      happens if the crystal's space group symmetry is lower than
!      its lattice symmetry.
!
RESOLUTION_SHELLS= 100 10 6 4 3 2 1.9
SPACE_GROUP_NUMBER=78
UNIT_CELL_CONSTANTS=57.39 57.39 106.9    90 90 90
OUTPUT_FILE=scaf8_all_merged.hkl
MERGE=TRUE FRIEDEL'S_LAW=FALSE
STRICT_ABSORPTION_CORRECTION=TRUE
INPUT_FILE= ../xds-1_2/XDS_ASCII.HKL
REIDX_ISET= -1  0  0  0  0  1  0  0  0  0 -1  0
INPUT_FILE= ../xds-2_1/XDS_ASCII.HKL
INPUT_FILE= ../xds-3_1/XDS_ASCII.HKL
INPUT_FILE= ../xds-1_4/XDS_ASCII.HKL
INPUT_FILE= *../xds-5_1/XDS_ASCII.HKL
```

Bold keyword/parameter pairs are required. Complete documentation at
xds.mr.mpg.de/html_doc/xscale_parameters.html

Output file **XSCALE.LP** shows level of systematic error **ISa** before/after scaling.


Example XDSCONV.INP

```
! UNIT_CELL_CONSTANTS= 10 20 30 90 90 90
! SPACE_GROUP_NUMBER= 96
! GENERATE_FRACTION_OF_TEST_REFLECTIONS=0.05

INPUT_FILE=XDS_ASCII.HKL
OUTPUT_FILE=temp.hkl CCP4_I+F ! or CCP4_I or CCP4_F or SHELX or CNS
FRIEDEL'S_LAW=FALSE ! store anomalous signal in output file even if weak
```

Bold keyword/parameter pairs are required. Complete documentation at xds.mr.mpg.de/html_doc/xdsconv_parameters.html

3 examples of single data sets (plots from *XDSGUI*)



XDSCC12: calculates $\Delta CC_{1/2,i} = CC_{1/2,with_i} - CC_{1/2,without_i}$

Three resolution ranges (blue=low green=medium red=high) - i refers to batches of width 1°

- find bad frame ranges
- radiation damage

xdsc12 for data from XSCALE

overall CC1/2: 99.784 nref= 29744 (but the overall CC1/2 is meaningless!)

<CC1/2>: 68.971 (frequency-weighted average of CC1/2 in resolution shells)

CC1/2 in resolution shells:

99.8	99.6	99.6	99.3	97.8	92.9	81.3	52.7	33.4	12.5
------	------	------	------	------	------	------	------	------	------

CC* in resolution shells:

99.9	99.9	99.9	99.8	99.4	98.1	94.7	83.1	70.7	47.1
------	------	------	------	------	------	------	------	------	------

frequency, i.e. number of unique reflections in resolution shells:

970	1741	2257	2669	3005	3336	3604	3904	4104	4154
-----	------	------	------	------	------	------	------	------	------

headings for lines starting with a,b,c:

a: <CC1/2> of each dataset:

a:	reflections of this dataset only	reflections of all datasets
----	----------------------------------	-----------------------------

a:	set	nref	with	without	delta	nref	with	without	delta
----	-----	------	------	---------	-------	------	------	---------	-------

b: delta-CC1/2 in resolution shells

c: # reflections for delta-CC1/2

a	1	27910	66.807	60.517	10.560	28672	66.486	60.349	10.249
---	---	-------	--------	--------	--------	-------	--------	--------	--------

b	30.906	29.045	30.431	27.944	34.663	27.727	23.134	14.772	13.087	5.266
---	--------	--------	--------	--------	--------	--------	--------	--------	--------	-------


c	912	1654	2173	2527	2803	3169	3394	3679	3849	3750
---	-----	------	------	------	------	------	------	------	------	------

a	2	25492	68.433	61.215	12.423	26174	70.905	63.097	14.130
---	---	-------	--------	--------	--------	-------	--------	--------	--------

b	27.445	18.627	5.799	14.503	25.465	32.858	31.704	23.753	13.638	-.254
---	--------	--------	-------	--------	--------	--------	--------	--------	--------	-------

c	858	1486	2064	2349	2517	2928	3153	3435	3317	3385
---	-----	------	------	------	------	------	------	------	------	------

output for two datasets
in P1 that scale quite
well



xdsc12 writes a sorted XSCALE.INP.rename_me

```
SPACE_GROUP_NUMBER=      1
UNIT_CELL_CONSTANTS=      53.00      59.93      61.34  99.975 105.633 119.747
OUTPUT_FILE= temp.ahkl
PRINT_CORRELATIONS= FALSE
SAVE_CORRECTION_IMAGES= FALSE
FRIEDEL'S_LAW= FALSE

! median of delta-cc1/2 "only" values=      11.492
! noise (MAD) of these values=      0.932
! median of delta-cc1/2 "all" values=      12.189
! noise (MAD) of these values=      1.941
! median of delta-cc1/2-ano "only" values=      -2.201
! noise (MAD) of these values=      2.699
! median of delta-cc1/2-ano "all" values=      -2.010
! noise (MAD) of these values=      2.614
! input files sorted by deltacc12_only (highest first):
! deltacc12 only / all:   12.4231   14.1301 deltacc12-ano only /all:   -4.9002   -4.6243
INPUT_FILE=../xds_14_2/XDS_ASCII.HKL
! deltacc12 only / all:   10.5600   10.2488 deltacc12-ano only /all:    0.4976    0.6043
INPUT_FILE=../xds_14_3/XDS_ASCII.HKL
```

output for two datasets
in P1 that scale quite
well

For important data, do not rely on automatic data processing

There are no perfect criteria for the quality of a data set - “Table 1” does not tell the whole story. Synchrotrons typically run multiple pipelines, and the user has to choose ...

What can go wrong in automatic data processing?

- does not handle radiation damage (we need to discard frames towards the end of the data set, but where should we cut?)
- does not handle shadowed areas of detector
- does not handle indexing problems (multiple lattices, ice, ...) flexibly
- does not optimize processing

Automatic processing with GlobalPhasing's *autoPROC*, CCP4's *xia2*, HZB's *XDSAPP* is rather reliable for good data! Difficult data sets typically benefit from human insight.

Manual processing with XDSGUI

- problems in phasing and refinement
often due to bad / wrong data processing
- visually inspect frames; mask shadows
- optimize parameters, frame range,
resolution cutoff ..
- presents tables as plots
- interfaces to XDS through its files, e.g. XDS.INP
- user – extensible / modifiable commands

Manual checks are easily performed with *XDSGUI*

- before running *XDS*: check frames (*adxv*, *Albula*, *XDSGUI* ...)
- **Frame** tab: FRAME.cbf for match obs/calc, and masked areas
- **IDXREF** tab: for ice rings / indexed vs un-indexed reflections
- look at **indexing in reciprocal space**:
 tools tab/Further analyses/spot2pdb-coot command
 see XDSwiki:XDSGUI#tools
- **statistics** tab: *XDSCC12* and *XDSSTAT*

Optimize data processing

- XDSwiki:Optimisation#Re-INTEGRATEing_with_the_correct_spacegroup.2C_refined_geometry_and_fine-slicing_of_profiles

XDSwiki:Optimisation#using_the_refined_values_for_beam_divergence_and_d_mosaicity_for_re-integration

- **tools** tab:
 - * “Saving and comparing good results” and “Optimizing data quality”. After changing parameters, run “JOB=DEFPIX INTEGRATE CORRECT”, compare and save if better/restore old if worse.
 - * “Further analysis”: Inspect indexed/unindexed spots in reciprocal space
- consider use of *StarAniso* if anisotropy is strong (i.e. visible)

... may make the difference between structure solved or not, interpretable or non-interpretable map, good or bad refinement, ...

Summary

Data processing is the crucial link between your experiment and your structure

- * garbage in – garbage out!
- * for important data, do not rely on automatic data processing alone
- * manual checks are easily performed with *XDSGUI*
- * try to optimize data processing – this converts noise to signal, and may enable structure solution, and/or improve refinement