

DLS-CCP4 Data Collection and Structure Solution Workshop  
29 October 2024

# MX at Diamond

Marco Mazzorana



# The UK national synchrotron

2003 Established

2007 First User

(7 beamlines of which 3 MX)



3 rounds of funding



Medical  
Research  
Council



Science and  
Technology  
Facilities Council



Research Complex  
at Harwell



The Rosalind  
Franklin Institute



CCP4



# Many different instruments...

33 beamlines

8 Science Groups

2 national cryo-EM facilities


12k+ users

100+ companies


Structural biology labs

MPL, crystallization facility XChem


UK-XFEL Hub



...to investigate the structure of biological systems



# The power of combining techniques



Super resolution  
cryo-light  
microscopy

Cryo X-ray  
tomography

Cryo-electron  
tomography

Atomic models  
from MX fit to  
averaged particles  
of two  
intermediates


Understand interactions and action across a large range of spatial and time domains (often weak and transient)

Reveal inner workings of the cell with a full dynamic picture of living organisms

Molecular description of cellular organization can transform fundamental science and understanding of disease processes



# What crystallography tells us about macromolecules



# Challenges and synergies

- PADS
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus



- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long Wavelength



- VMXi
- I24
- I03
- Screening
- Data collection
- CLIII

*In situ*



- Long wavelength
- Multi-axis strategies
- Inverse beam
- Wedged MAD
- Pipelines

Experimental Phasing



- *In situ*
- Containment Level II & III
- Schedule 5

Biological Containment




- Humidity control
- Spectroscopy
- Tomography
- Dynamics
- Fluorescence

Complementary information



- Diffraction image processing for SR
- Electron and XFEL data
- Collaborative effort
- Led and managed by DLS

Dials Project




- Experiment database
- Synchweb
- Data archive
- Data access
- Data reprocess

Data Management




# Dealing with the smallest crystals: microfocus


3 microfocus beamlines



I04




I24



VMXm

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus







# Microfocus for time-resolved studies

Sample delivery for serial crystallography (SSX and SFX)


Liquid jets/ LCP extruders



Acoustic droplet ejection




Fixed target



Silicon nitride chips  
25,600 positions/chip


- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




# Microfocus for time-resolved studies

Fixed-target at I24




Bespoke mount in I24  
Fast stages



Full chip rastered in 6-20 mins  
10-50 ms/frame (+ pause at each frame)  
Typically ~25 chips/12 hours


Multi-crystal structures (from stills)



1.5 Å structure  
from 9000 crystals

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




# The UK-XFEL Hub



- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




tools and strategies to enable dynamic structural biology facilitating:

- (time-resolved) serial structural biology experiments via sample preparation, delivery, data collection, and processing
- the transfer of methods between XFEL, synchrotron, and/or cryo-EM sources
- access to, and data collection from complementary facilities (DLS, CLF, eBIC)






# Extending the time range





- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus




# The Diamond-II project

4th generation synchrotrons and XFELs have hotter beams



# Smaller than visible

## VMXm: the first nano-focus beamline




< 500 nm beamsize

Integrated SEM


In vacuum

Tunable 7-25 keV (2.0-0.57Å)



- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*

Microfocus








# Smaller than visible


VMXm: the first nano-focus beamline

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*


Microfocus



OAV




SEM



extremely low background

# Smaller than visible

VMXm: the first nano-focus beamline




27 crystals

1.5-2.0 deg/data per crystal

Rigid body+restrained refinement

Dials, CCP4, Refmac, Coot




Summary of merging statistics

	Overall	Low resolution	High resolution
Resolution (Å)	36.57 - 1.50	36.58 - 4.07	1.53 - 1.50
Observations	116543	7779	1406
Unique reflections	25790	1405	833
Multiplicity	4.5	5.5	1.7
Completeness	87.23%	90.30%	57.21%
Mean I/ $\sigma$ (I)	5.1	13.4	0.7
Rmerge	0.195	0.128	0.775
Rmeas	0.217	0.140	1.017
Rpim	0.088	0.053	0.646
CC1/2	0.969	0.970	0.319

- Tuneable
- Rapid
- Variable focus
- High flux
- *In situ*


Microfocus





# Not just photons – The HeXI beamline

## High energy electrons for structure determination

- HeXI will target 300 nm to 3  $\mu\text{m}$  crystals with high energy electrons
- Bridge the gap between Commercial electron Microscopes and MX beamlines




## Project Milestones




- Funded by the Wellcome ELS grant
- Conceptual design approved
- Technical design submitted
- Electron source tender submitted
- Commencing construction

**First Light by 2026**

Operation through the  
"Dark Period"



## Dedicated MeV electron diffractometer



New Grant will Electrify Life Sciences  
with Breakthrough Technologies






# The need for longer wavelength

Biologically relevant elements

	e (keV)	$\lambda(\text{\AA})$
P	2.146	5.779
S	2.472	5.016
Cl	2.822	4.393
K	3.607	3.437
Ca	4.038	3.070

Xray absorption edges outside the capabilities of standard beamlines





- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long  
Wavelength



# Challenges and solutions



- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long  
Wavelength




Long wavelength increases the angle of diffraction  
Need a curved detector to measure high resolution spots


# Challenges and solutions

Absorption by air of  
long wavelength X-rays

In vacuum beamline



Bespoke sample holders




- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long  
Wavelength







# Results: metal identification



Peptidyl transferase center



Decoding center




- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems


Long  
Wavelength




# Extending the element 'palette'




ribosome 70S




Na/K transporter




FAD-binding



SLC7 transporter



PETase



- Tuneable beamlines
- I23 dedicated
- Smart data collection strategies
- Expert systems

Long  
Wavelength





# High throughput



BART sample changer:  
37 unipucks  
597 samples  
Exchange rate approx. 20 sec/sample

## EIGER2 XE 16M

Pixel size:  $75 \times 75 \mu\text{m}^2$

18,093,576 pixel

560 frames/sec (3600 frames in <7 sec)



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput






# Logistics



## Weland/Randex automated store

high capacity (up to 272 dewars)  
barcode tracking (dewars and pucks)  
logistics integrated in Synchweb/ISPyB



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



# Responsive/dynamic scheduling

2019

2020

[illegible]

Run	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100	
Run 2	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 3	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 4	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 5	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 6	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 7	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 8	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 9	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 10	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 11	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 12	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 13	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 14	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 15	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 16	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 17	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98	99	100
Run 18	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59																																									

[illegible][illegible][illegible]


- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



# Unattended data collection

Optimized for speed




- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput




Xray-centering and strategies





# Automated processing pipelines



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput





# HT fragment screening

## Conventional screening

Preselection of best compounds

### Conventional screening

- Searching for potent molecules
- Complex molecules → Low probability





## Fragment-based approach

Screen by crystal

### Fragment screening

- Guaranteed binding – but weak
- Potency through chemical elaboration



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)


High Throughput




# HT fragment screening

## XChem Lab at Diamond


Fragment libraries




## Beamline I04-1



~700 datasets/day, unattended



Protein crystals in plates



- PADs
- Automation
- Pipelines
- Serial Crystallography
- Access (rapid/freq)

High Throughput



1 crystal – 1 fragment

Soak/harvest/collect data  
(up to 1000 crystals/week)

Routine since 2016

March 2020 COVID-19 targets

# Multiple access ways

## Academic

Peer-review access based on quality of proposal

- Free access
- Publishable experiments

Access modes:

- **BAG** (2 year programme)
- **Single application** (6 months)
- **Rapid access** (apply anytime for one-off experiments)

MX beamlines, XChem



## Proprietary (industry and academia)

Organized via Industrial Liaison Office

- No need to publish
- Allocated slots throughout the schedule
- 10% of allocated beamtime

Services provided:

- Direct and remote access
- Unattended data collection
- Mail-in data collection and analysis service
- XChem fragment screening



 [instruct-eric.org](https://instruct-eric.org)  [admin@instruct-eric.org](mailto:admin@instruct-eric.org)  [@instructhub](https://twitter.com/instructhub)

 **Instruct-UK- Harwell Campus**






# Access to MX

## Cryo



standard SPINE pins  
I23 specific holders  
VMXm – EM grids



20 µm

## SSX (I24, VMXi and XFEL Hub)

Fixed targets  
Thin films  
Extruders  
Acoustic drop ejection  
Tape drives (in dev)



## Room Temperature



VMXi dedicated beamline  
I03 - biological containment

# Crystallisation facility


>£1.4M invested in last 5 years (DLS/ RCaH)  
Discussion over future partnership with ISIS  
Rockmaker and Rockmaker web (Formulatrix)

## Future Plans

- Replacement of 4 °C imager
- Imager for nano crystals- VMXm/ microED
- Upgrade chemical storage including Echo dispensing
- Long-term technical support (1.0 FTE) for facility



Total number of plate throughput  
by Financial Year (Apr-Mar)



**2018**  
Formulatrix R1000  
LCP plate Capacity

**2018**  
Scorpion

**2019**  
VMXi  
Workflow

**2022**  
Formulatrix R1000  
LCP + MFI optics

**2022**  
Formulator  
+ LCP Mosquito

**2023**  
Mass  
Spectrometer

**2024/25**  
Imaging for nano xtals  
tracking system



[instruct-eric.org](https://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)



# *In situ* data collection

Fully automated handling of *in-situ* plates at 20 °C for data collection (100s of datasets/hr)

~10x10  $\mu\text{m}^2$  pink beam  $3 \times 10^{13}$  ph/s, 16 keV

Machine learning crystal location and selection


Users can remotely queue plates for unattended data collection

Optimisation of crystallisation conditions as well as high quality room temperature structures


Multi-crystal data merged using **dials-multiplex** and presented in ISPyB with downstream analysis




*Crystal hotel to the beamline using robotic arm*



*Up to 60° rotation of plates*



*Machine learning crystal identification & marking for data collection*



*Automated multicrystal data processing*

**Plate arrival to processed data can be within a few hours**



[instruct-eric.org](http://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)





# Membrane protein laboratory (MPL)

Contact: [mpplab@diamond.ac.uk](mailto:mpplab@diamond.ac.uk)

Andrew



James



Harish



Claudia



Peter



Rebecca



Charlie



- Funded by Wellcome - free at point of access
- Access: **Rapid proposal** & **Instruct-ERIC**
- **Platforms:**
  - HTP cloning and expression screening
  - Purification Optimisation
  - Sample preparation for cryo-EM, crystallisation and SAXS
  - Biophysical Characterisation
  - Data collection & structure solution



Research Complex  
at Harwell



[instruct-eric.org](https://instruct-eric.org)




[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)




# eBIC access through INSTRUCT



CryoEM SPA

Established Since start




Cryo-Electron Tomography

Established & Developing




Chameleon & Electron Diffraction

Developing




Leica CryoCLEM


Developing




Nayak *et al.*, 2023 NSMB



McLaren *et al.*, 2023 Nat.Microbiology



Radiation hard 4KHz readout



Humidity control chamber



[instruct-eric.org](https://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)



instruct  
ERIC





diamond

# RELION 4 Pipeline

with new web interface developed outside of SynchWeb Introducing PATo

(Dan Hatton, Anna Horstmann, Stephen Rigg, Guilherme De Freitas)




**Data storage**  
18 TB/SPA visit

Integrated data transfer and triggering the pipeline - the local contact needs to input the dose per frame and everything else is automated (relion, cryolo, ice breaker) And it is Live!



# B24 correlative imaging



ation

Computerised Tomography  
1-2 mm

Fluorescence Microscopy  
250 nm

Cryo-Soft X-ray Tomography  
25 nm

Cell Electron Tomography  
5 nm

Single Particle Electron Microscopy  
0.2 nm

X-ray Crystallography  
< 0.1 nm



- Cryo soft X-ray (cryoSXT)
- 3D X-ray diffraction imaging




# Molecular Biophysics



2x ITC  
2x AUC  
3x RT-PCR  
Biacore T200  
SEC-MALS

~£1M invested in last 5 years (RCaH/ DLS)

**Future Plans**  
Echo dispensing



**2018**  
Stopped Flow



**2019**  
NanoDSF



**2020**  
Clariostar



**2022**  
Fida1



**2023**  
Mass Spectrometer



**2025**  
Echo dispensing



[instruct-eric.org](https://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)

# B21 – bioSAXS/WAXS

Source: bending magnet

Flux:  $4 \times 10^{12}$  ph/s

Wavelength:  $0.9464 \text{ \AA}$  (13.1 keV)

Sample to SAXS detector distance: 3684.6 mm

SAXS Detector: EigerX 4M (Dectris)

WAXS Detector: EigerX 1M (Dectris)



SAXS Q range:  $0.0045\text{--}0.34 \text{ \AA}^{-1}$

**NEW** WAXS Q range:  $0.8\text{--}3.4 \text{ \AA}^{-1}$

## Upgrade plans:

New Agilent HPLC online by December 2024 with  
Integrated GDA/HPLC controls in 2025

Robot Arm for solid and viscous samples installed by late 2025



[instruct-eric.org](https://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](#)




instruct  
ERIC






diamond



# INSTRUCT internships – a success case



**129 fragments bound to the protein**



[instruct-eric.org](http://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)




# INSTRUCT internships – a success case

## AI-based automated design

Syndirella  
HIPPO

590 hypothetical  
elaborations (8000€)




## Robotic synthesis

CAR

454  
elaborations


252  
products



## Biophysical assays


Creoptix

6 pure compounds with  
 $\mu\text{M}$  affinity for the protein



## High-throughput crystallography

16 crystal structures  
from 3 scaffolds  
(crude reaction mixtures)



[instruct-eric.org](http://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)



# Instruct-UK at Harwell Campus



**Diamond Light Source**

- **Fragment Screening**
- Electron microscopy
- X-ray diff. & Bio-SAXS




**Research Complex at Harwell**

- Molecular biophysics
- *Membrane protein expression*
- Membrane protein crystallisation
  - Crystallisation



Harwell Campus

**Admin Contact:**  
Jodie Lavender




**Sci/Tech Contact:**  
Martin Walsh



**Scientific Contact:**  
Gwyndaf Evans



**Technical Contact:**  
Gemma Harris




**Scientific Contact:**  
Andrew Quigley



**Scientific Contact:**  
Frank von Delft



**Technical Contact:**  
Halina Mikolajek



**Scientific Contact:**  
Karen Davies



**Technical Contact:**  
Daren Fearon



[instruct-eric.org](https://instruct-eric.org)



[admin@instruct-eric.org](mailto:admin@instruct-eric.org)



[@instructhub](https://twitter.com/instructhub)





# The team



dls.mx



mx-usersupport@diamond.ac.uk

