





Radiation Damage


CCP4/DLS MX Workshop:




Virtual lecture
18 November 2024



elspeth.garman@bioch.ox.ac.uk




Radiation Damage The Plan:



- **Why cool? Radiation damage**
- What are the symptoms?
- What is it?
- Why do we care? Effect on MAD/SAD.
- How do we estimate the Dose?
- The limits
- What can you do to minimise it?

I24, Diamond, *in situ* data collection from a
Bovine Enterovirus 2 crystal, room temperature, 0.5 s
20 μm x 20 μm beam



Axford *et al.*,
Acta Cryst D (2012) 592





Beamline logo I24, 5 μm beam, large crystal!
(Gwyndaf Evans *et al.*)

ROOM TEMPERATURE



Intensity decrease
with dose
especially at high
Bragg angles


Loss of
diffraction

Incomplete data
from crystals at
RT




Radiation damage

Primary 
Secondary 



Protein: direct
Solvent: indirect





PRIMARY; inevitable, a fact of physics! Neutralise it?



SECONDARY, can we control it?

Radiation damage





Significantly reduced at
100 K: time factor of ~70

[Nave and Garman *JSR* (2005), **12**, 257-260].

SP445: Data Quality


[T.Schneider]



PRIMARY; inevitable, a fact of physics! Proportions?



SECONDARY, can we control it?



Radiation Damage

The Plan:

- Why cool? Radiation damage
- **What are the symptoms?**
- What is it?
- Why do we care? Effect on MAD/SAD.
- How do we estimate the Dose?
- The limits
- What can you do to minimise it?




ROOM TEMPERATURE

Intensity
decrease
with dose
especially at high
Bragg angles


Loss of
diffraction

Incomplete data
from crystals at
RT

$$2d \sin \theta = n \lambda$$





Intensity decay at RT and 100 K




Radiation damage and dose limits in serial synchrotron crystallography at cryo- and room temperatures. De la Mora, Coquelle *et al.* (2020) *PNAS*



RadDam signatures in reciprocal space 100 K




(a)



(b)




(c)




(d)

RadDam signatures (100 K MX)



A



B




C




SSRL, 19 hours, 9.1, 1998

1995 onwards: 100 K
BUT THEN, 1999:



2002:
30 μ m
beam




[Tassos Perakis]



Also observe
spectral changes

Iron containing protein, ESRF

Garman & Owen (2006), Acta D62




At 100 K

Intensity
decrease

Loss of
diffraction



Incomplete data
from crystals



Happens during 1 dataset at 100K for many crystals

Unit cell volume expansion

Increase in Wilson B factors, Rmerge

Increase in mosaicity

'GLOBAL' damage

ESRF 2000:

1×10^{12} ph s⁻¹ into
100μm square slits

Australian synch.

3×10^{13} ph s⁻¹ into

50μm × 70μm

[10^{14} ph/s/ 100μm²]

Diamond Light Source


(I24): 3×10^{12} ph s⁻¹

into 7μm × 6μm

[7×10^{14} ph/s/ 100μm²]



RadDam signatures in reciprocal space 100 K



Raimond
Ravelli
1968-2023


What global damage metric should we use and against what should we plot it?

- I_n/I_1

- Not $I/\sigma(I)$

- Scaling B-factors?

- An R_{meas} type measure?



[Diederichs 2006 Acta D59, 903]

To monitor the damage rate, compare results from different experiments and avoid exceeding dose limits, we can't use time or image number.

We need a **DOSE** estimator

$$\text{Dose} = \frac{\text{energy absorbed}}{\text{unit mass}} = \frac{\text{J}}{\text{kg}} = \text{Gy (gray)}$$

Fundamental metric against which to plot damage indicators.

Cannot be measured, can only be estimated.

Takes care of the physics but NOT the chemistry.

RADDOSE-3D: www.raddo.se

RADDOSE-3D

V1: Zeldin, Gerstel, Garman (2013) *J.Appl.Cryst*,

V2: Bury et al (2018) *Protein Science*

V3: Dickerson and Garman (2021) *Protein Science*


V4: Dickerson et al. (2024) *Protein Science*

Intensity Decay at 100K


Normalised Intensity vs Dose:

apoferritin

[Owen et al 2006, PNAS]



[Diederichs 2006 Acta D59, 903]



R_d : pair wise R factor between identical and

symmetry-related reflections occurring on different diffraction images


Coefficient of sensitivity α

change in relative isotropic B factor:

$$s_{AD} = \Delta B_{rel} / 8\pi^2 \Delta D$$

(e.g. HEWL@100 K = 0.012 Å²/Gy)


[Kmetko et al 2006, Acta D62, 1030]



Real Space: Specific Damage

CEWL
lysozyme
Fo₄-Fo₁

Raimond
Ravelli





Real Space: Specific Damage

- Specific damage occurs in a predictable hierarchy in proteins
- Reduction of metallocentres
- Breakage of disulphide bonds
- Asp and Glu decarboxylation
- Cleavage of S—C bond in MET
- Rupture of covalent bonds to heavier atoms:

C-Br, C-I, S-Hg

Note that if this were all due to primary damage alone, damage would be in order of absorption cross sections of atoms, which it is not.



Difference map $Fo_4 - Fo_1$


Weik et al. 2000, Burmeister 2000
Ravelli & McSweeney 2000
Bury et al., JSR 2015

CO₂ from decarboxylation

Raddam induced decarboxylation of Asp128.

0.5:0.5 occupancy of C and 2 O atoms involved


N9 influenza neuraminidase, PDB: 6HFC




Resolution: 1.29 Å

R-Value Free: 0.164

R-Value Work: 0.141




ALSO:
Atomic B-factors increase:




Overall summary of RD effects

Global radiation damage


Decreased reflection intensities




Unit cell expansion



Scaling B factors




Increased mosaicity




Specific radiation damage

Chemical changes


Reduction of metal ions




Reduction of disulphide bonds




Glutamate / aspartate decarboxylation




Side / main chain motion





Why we cool...?


The Plan:



- Why cool: radiation damage
- What are the symptoms?
- **What is it?**
- Why do we care? Effect on MAD/SAD.
- How do we estimate the Dose?
- The limits
- What can you do to minimise it?

PHYSICS of the interaction of X-rays with crystals.

A) Diffraction






B) Absorption = Energy loss



N.B. $> 90\%$ of the beam does not interact at all,
but goes straight through.

Interaction of X-rays with biological samples


At $E_x = 12.4$ keV, 100 μm protein crystal, only 2% of beam interacts.



Photoelectric Cross Sections (barns/atom) at 13.1 keV

[1 barn= 10^{-28}m^2]


A few heavy atoms can
make a big difference.



Photoelectric Cross Sections (barns/atom) at 13.1 keV

[1 barn = 10^{-28}m^2]


A few heavy atoms can
make a big difference.




Photoelectric Cross Sections (barns/atom) at 13.1 keV

[1 barn=10⁻²⁸m²]

A few heavy atoms can
make a big difference.



Beam absorption ($\lambda=0.95\text{\AA}$, 13 keV) by a Zn containing 30 μm protein (metallo- β -lactamase) crystal




Damage: the Radiation Chemistry

1) INDIRECT RADIATION DAMAGE :

radiolysis of water

at room temperature
and neutral pH:




OH thought not to be
mobile in glasses
below 110K

(Owen et al Acta D
2012)



Hiroki, A. Pimblott, S. M.
LaVerne, J. A. (2002)
J Phys Chem A **106**,
9352-9358

Radiation Chemistry in action: Nitrate scavenger



Why we cool...?

The Plan:




- Why cool: radiation damage
- What are the symptoms?
- What is it?
- **Why do we care? Effect on MAD/SAD.**
- How do we estimate the Dose?
- The limits
- What can you do to minimise it?

Why do we care?

A) We don't get all the data we need!

B) Effect on MAD/SAD phasing methods.




- Failure of structure determination

(Multi-wavelength anomalous dispersion MAD, SAD)

due to creeping non-isomorphism –

- a) cell expansion and
- b) movement of molecule in unit cell
- c) structural changes DURING experiment.

i.e. **MAD/SAD phasing** signals (<5%) washed out completely




e.g. for a 0.5% change in all 3 unit cell dimensions of 100Å, reflection intensity changes by **15% at 3Å**

Why care?

C) Radiation damage can affect our biological results

- Metallo-proteins often photo-reduced during the experiment [e.g. PSII, Yano et al, PNAS (2005)]
- Decarboxylation of Glu is sometimes part of the protein mechanism, but is indistinguishable from radiation damage at the synchrotron
- X-ray induced structural changes were initially misleading in studies of intermediates e.g. Bacteriorhodopsin: orange species





Takeda et al, Crystal structure of the M intermediate of bacteriorhodopsin...JMB (2004),
Wickstrand, et al., Bacteriorhodopsin: Would the real structural intermediates please stand up?,
Biochim. Biophys. Acta (2014)

Why we cool...?

The Plan:

- Why cool: radiation damage
- What are the symptoms?
- What is it?
- Why do we care? Effect on MAD/SAD.
- **How do we estimate the Dose?**
- The limits
- What can you do to minimise it?



DOSE Limit Postulate

(Henderson 1990):

- There is a MAXIMUM dose for MX postulated as 20 MGy
(Energy absorbed/unit mass: Joules/kg = Gy)
tolerated by protein crystals: depends only on PHYSICS of situation.
Henderson: 20 MGy
Experimental: 43 MGy to $I_{0.5}$ – but don't go lower than $I_{0.7}$, 30 MGy*
- Crystal might not reach that limit due to chemical factors, but it is unlikely to last BEYOND the limit.
- Need to be able to calculate the DOSE conveniently: **RADDOSE-3D**

RADDOSE-3D

V1: Zeldin, Gerstel, Garman (2013) *J.Appl.Cryst*,

V2: Bury et al (2018) *Protein Science*

V3: Dickerson and Garman (2021) *Protein Science*


V4: Dickerson et al. (2024) *Protein Science*

*Owen, Rudino-Pinera, Garman (2006), PNAS


Make the dose calculation convenient for MX (include solvent contribution in mM and heavy atoms explicitly)

To find the energy deposited per unit mass in the crystal, need to characterise two things:

The beam




The crystal




Calculating Dose (*RADDOSE-3D*)


Crystal Characteristics



absorption coefficients
e.g. apoferritin: 0.406mm^{-1}
holoferritin: 1.133mm^{-1}



200 μm





Number of Amino Acids

‘HA’ atoms per monomer, e.g. S, Se, Hg

Solvent - concentrations of components, e.g. Na⁺, Cl⁻

Calculating Dose (*RADDOSE-3D*) Beam Characteristics




PIN diode to measure flux (ph/s)

Calculating Dose (*RADDOSE*)

[R-v1: Murray, Garman & Ravelli, JAC 2004
R-v2: Paithankar, Owen & Garman, JSR 2009,
R-v3: Paithankar & Garman, Acta D 2010]

Crystal Characteristics



Beam Characteristics



e.g. absorption coefficients, μ_{abs} , at 12.7 keV
apoferritin: 0.406mm^{-1}
holoferritin: 1.133mm^{-1}

RADDOSE-3D

- TIME- and SPACE-resolved modeling of dose distributions in MX in Java, replaces RADDOSE:
- Full 3-D simulation of dose absorption by the crystal
- Can deal with multiple wedges of data and different energy beams (e.g. MAD)
- Models beam as Top-Hat or Gaussian **or can use measured experimental profiles**
- Engineered for easy extendibility:
can use any crystal shape
- On server: www.raddo.se (!!)
- On github.com/GarmanGroup/RADDOSE-3D



[Zeldin, Gerstel, Garman *JAC* (2013)


Bury, Brooks-Bartlett, Walsh, Garman, *Protein Science* (2018)]

RADDOSE-3D New GUI

Download from: <https://github.com/GarmanGroup/RADDOSE-3D>

Versions for a PC and for Linux, no version for a Mac (anyone interested?!):


RD3D_GUI_Windows.exe and RD3D_GUI_linux



Josh
Dickerson

Dickerson *et al*
Protein Science
(2024) 33.e5005

Dose distribution vs exposure strategy with RADDOSE-3D



0.0001 MGy (grey),
5 MGy (green),
10 MGy (light blue),
20 MGy (dark blue),
30 MGy (red),
360° rotation

20 μm
translation

Helical,
0.2 $\mu\text{m}/^\circ$

30 μm
offset


30 μm offset
and helical
0.2 $\mu\text{m}/^\circ$

Helical strategy improvement:
Flot *et al* (2010) *JSR* **17**, 107

100 μm , 200 μm , 100 μm .

Gaussian beam (FWHM: 20 μm x 20 μm), 12.4 keV, 5×10^{11} ph/s,
1 mm x 1 mm rectangular collimation: full crystal bathed in beam


Zeldin *et al.* *JSR* (2013), Bury *et al.* *Protein Science* (2018)



Radiation Damage

The Plan:

- Why cool: radiation damage
- What are the symptoms?
- What is it?
- Why do we care? Effect on MAD/SAD.
- How do we estimate the Dose?
- **The limits**
- What can you do to minimise it?



Typical MX experiment

1 MGy/s absorbed by a 100 μm cubed metal free crystal in a $100 \times 100 \mu\text{m}^2$ 12.4 keV (1 Å) X-ray beam of flux: 10^{13} photons s^{-1}

MX at RT: ~0.4 MGy experimental dose 'limit' (dose to $0.5 I_n/I_1$) reached in ~0.4 s

MX at 100 K: 30 MGy experimental dose 'limit' (dose to $0.7 I_n/I_1$) reached in ~30 s:
4th generations sources $\ll 1\text{s}$,

MX at XFELs: damage before destruction ?
<80 fs




2-10 Gy, dies
in 10-30 days
Radiotherapy on
human brain
tumour max 60 Gy

RT



40 consecutive frames per position

Equivalent-dose datasets
(19818 positions at each dose)



De la Mora, Coquelle, .Holton...Garman, *et al.* Weik
(2020) PNAS + Suppl Material


RadDam signatures in reciprocal space (RT MX)



J Crystal Growth 90 (1988) 259-272

Radiation damage and dose limits in serial synchrotron crystallography at cryo- and room temperatures. De la Mora, Coquelle *et al.* (2020) PNAS


- At 100 K, specific damage effects onset before global damage
- Both classes of damage common in protein crystallography (MX)



[1. Henderson Phil Trans RS B 1990

2. Owen *et al* PNAS 2006]



Time to dose limits at current synchrotrons



Average diffraction weighted doses (DWD) were calculated for a $(50 \mu\text{m})^3$ lysozyme crystal (grown in 100 mM NaAc and 1M NaCl, solvent fraction 38%) rotated 360° in the flux within the FWHMs of the various beams.

Why we cool...?

The Plan:




- Why cool: radiation damage
- What are the symptoms?
- What is it?
- Why do we care? Effect on MAD/SAD.
- How do we estimate the Dose?
- The limits
- **What can you do to minimise it?**



1) Sacrificial crystal to characterise damage rate



Diffraction fading



I/I_1



Dose = absorbed energy (J) / mass (kg)


2) During the experiment:

- Very hard to see spot fading by eye on Eiger detector as all reflections are partials. Crystal usually diffracts further out than you can see.
- Eiger has zero dark current background.




RadDam signatures in reciprocal space


Undamaged




Damaged



Scale and R_{merge} vs batch



Scale and R_{merge} vs batch



DLS

Ed Lowe

3) DOSE AWARE data collection

Some beamlines worldwide are running RADDOSE-3D
'underneath' to allow collection by DOSE, not TIME

e.g. RADDOSE-3D is now embedded within Blu-Ice at SSRL so runs on all the MX
beamlines there.

Also on I04 at DLS

(David Aragao)

The screenshot displays the Blu-Ice software interface for BL14-1. The main window shows a diffraction pattern with a central bright spot. Overlaid on the left is a yellow-bordered box containing the 'Run Dose Estimate (MGy)' panel, which displays the following information:

- Room Temp Limit: 0.38
- Predicted: 1.81
- Exposed: 0.00

The right side of the interface contains several control panels:

- WebSite Strategy:** Collect, Pause, Abort, Advx Autoload.
- Current Position:** Phi: 18.71 deg, Distance: 500.00 mm.
- Run Dose Estimate (MGy):** Room Temp Limit: 0.38, Predicted: 1.81, Exposed: 0.00.
- Optimize Beam:** A list of test parameters (test_1_00001 to test_1_00024) with values for Phi and Distance.
- Run 1 (inactive):** Prefix: test, Directory: data/scottm, Beam Size: 130.0 x 80.0 um, Detector: normal, Distance: 500.000 mm, Beam Stop: 43.488 mm, Axis: PHI, Delta: 1.00 deg, Transmission: 100.000 %, Time: 1.0000 s.
- Frame:** Start: 1, End: 1000, Phi: 18.71 deg, Inverse Beam: 999999.00 deg, Wedge: 8543.020 eV, Energy: 8543.020 eV.
- Run Time:** 00:16:42.
- Pause Score Update:** A checkbox.
- File Name:** /data/debanu/ab2022/ab_xtx/20220705/collect/xtx11_31/xtx11.
- Scripts reports:** A log of script execution.
- System Status:** Cryostat temperature 202 K, Spear Current: 499.803 mA, Abort, User: Active, Shutter: closed, 12:25:22 PM.

Dose and Exposure

Number of Images

3600

☒ Set Target Dose

☐ Set Target Exposure

Exposure Time

0.0020

s

Total Exposure Time

7.2

s

Dose / Dataset

0.5

MGy

First Image Number

1

Rules of thumb at 100 K (from my colleague Ed Lowe with much experience):
Spread 10 to 15 MGy over 360° to 720°
Max dose: 7 MGy if have 1.4 Å diffraction, can allow more dose for lower resolution
Diminishing returns beyond 15 MGy
AVOID TEMPTATION TO INCREASE THE EXPOSURE TIME or DOSE




4) BACK SOAK to REMOVE non-specifically bound heavy atoms

A few heavy atoms in the solvent can make a BIG difference to the absorption cross section and this the dose rate for the SAME flux.




e.g. Cacodylate buffer (arsenic, mass 75 cf selenium = 79) [Tip: avoid cacodylate!]

e.g. a brominated DNA-protein complex will radiation damage much faster than a native crystal and will debrominate during data collection [Ennifar et al, Acta Cryst D (2002) 1263-1268].

5) If possible, use a top-hat like beam



Imaged beam
I24, DLS,
Axford, Owen



Imaged beam
P14, PETRA III,
Bourenkov,
Schneider


Garman & Weik COSB (2023)

6) For data collection: consider optimising the incident energy

Monte Carlo simulations: CdTe detector, peak in DE@ 26 keV
diffraction efficiency (DE) = diffraction/dose


Photoelectron escape

Experimental verification: Storm *et al* IUCrJ 2021



7) Exclude damaged images

High resolution data:




555 images

(100 K, 1 Å resolution)

high, medium, low, high2

must exclude images
showing radiation
damage from the dataset

588 images




8) During refinement:

There are <30 ‘radiation damage series’ in the PDB.

Can we give an isolated deposited PDB file a **YES!!** ‘radiation damage index’?




Kathryn Shelley




The B_{Damage} metric

- B_{Damage} is B -factor corrected for packing density



RABDAM (now in the CCP4 computing suite)

- RABDAM calculates B_{Damage} for all selected atoms in any standard format PDB file
- B_{Damage} highlights expected sites of specific radiation damage
- RABDAM provides several useful outputs to aid radiation damage evaluation




RABDAM (now in the CCP4 computing suite)

Get a 'radiation damage index, B_{net} ' for a 100 K


Use CCP4 program RABDAM. Deconvolutes packing density from B -factor to give B_{Damage}

Uses B_{Damage} to compute B_{net}




The B_{net} metric

- B_{net} is calculated from the distribution of the B_{Damage} values of Asp $\text{O}\delta$ and Glu $\text{O}\epsilon$ atoms
- Equal to the ratio of the area under the curve either side of the median B_{Damage} value of all protein atoms




$$B_{\text{net}} = \frac{\text{Area RHS}}{\text{Area LHS}}$$



The B_{net} metric


- B_{Damage} is a per-atom metric
- The B_{net} metric is a derivative of B_{Damage} that summarises the total extent of specific radiation damage suffered by a PX structure in a single value




Analysis of the Protein Data Bank: *B*_{net}

- Most of the structures with the highest B_{net} values show clear signs of radiation damage in their electron density. Above $B_{\text{net}}=3$ look at structure before you use it!



3IWU ($B_{\text{net}} = 78.2$)





5G28 ($B_{\text{net}} = 66.3$)




5WUC





1V70




5FXL





6Q5R





5XQP





3A07





3S8S





6BKL



3UX1




2XMK



FURTHER READING:

- ‘Beginner’s guide to Radiation Damage’ Holton, (2009)
JSR **16**,133-142
- General summary in: Garman, *Acta D* (2010) **66**, 339-355
- Garman and Weik, Chapter 20 in ‘*Protein Crystallography: Methods in Molecular Biology*’ (2017) **1607**, 477-489



14th March 2022

Quantifying and comparing radiation damage in
Protein Data Bank

Kathryn L. Shelley ^{1,2} & Elspeth F. Garman ¹

- Shelley and Garman, *Acta D* (2024) **80**, 314–332

The Crystallographer's DILEMMA:



Rate of damage
versus diffraction
intensity

This was a lot of material delivered quickly so...



**Please consider
publishing in
IUCr Journals
Acta D or F:
see next slides!**



**Questions very welcome, now and
when we finally meet**

elspeth.garman@bioch.ox.ac.uk

EFG by Lorna
Dougan,
Leeds,
LOOminaries



Structural Biology


Acta Crystallographica Section D



journals.iucr.org/d/



@ActaCrystD



Society journal established in 1993

Covers all aspects of structural biology

New structures of biological importance from any technique

Reports on the latest software, equipment and methods

International Editorial Board and trusted expert reviews

Easy submission with technical editing upon acceptance

Free to publish – no submission fees or page charges

Open access option – free if your institution is part of a deal



Structural Biology Communications


Acta Crystallographica Section F



journals.iucr.org/f/



@ActaCrystF



Society journal established in 2005

Covers all aspects of structural biology - concisely and rapidly

Biological structures determined by any technique

Brief reports on the latest software, equipment and methods

International Editorial Board with trusted expert reviewers

Early career researchers are encouraged and highlighted through interviews with our Editors

Free to publish - no submission fees or page charges

Open access option - free if your institution is part of a deal

Never too young to
start learning
crystallography...

Madeleine, aged 9 months, my
granddaughter reading BCA
Crystallography News!


Crystallography News

British Crystallographic Association



Issue No. 162 September 2022





ISSN 1467-2790



AlphaFold's implications for structure determination

AlphaFold: the Start of Something Bigger?	p6	CryoEM Book Review	p20
Spring Meeting 2023	p12	ECA Max Perutz Prize	p21
National Crystallography Service	p14	UKRI Infrastructure Fund	p22
Meeting Reports	p16	News from the CCDC	p24

RD2: Dec 2001 RD3: Nov 2003 RD4: March 2006 RD5: March 2008



JSR, Nov 2002 (8)

JSR, May 2005 (9)

JSR, Jan 2007 (14)




JSR, March 2009 (8)

RD6: Mar 2010

RD7: Mar 2012

RD8: Apr 2014

RD9: Mar 2016



JSR, May 2011 (10)

JSR, Jan 2013 (6)

JSR, March 2015 (8)


JSR Jan 2017 (8)

RD10: Sep 2018



JSR, July 2019 (9)

RD11: Sep 2021



JSR, July 2021 (6)

RD2 to 11: 86 papers published in JSR Special Issues to date

STOP PRESS!! RD12 3-5th June 2025, PSI, Switzerland

See: <https://indico.psi.ch/e/rd12>

RADDOSE-3D

TEST our new GUI!!

To run RADDOSE-3D for MX, SMX or SAXS (which ever you like!)

Step 1: Download and unzip the RADDOSE-3D GUI from:

<https://github.com/GarmanGroup/RADDOSE-3D>

There are versions for a PC (Windows_release.zip) and for Linux (Linux_release.zip).

If you have a MAC, there is no new GUI yet, but you can run a limited capability RADDOSE-3D from the WWW site:

raddo.se


(click on ‘manual interface’ and run the test example first. Then edit the input for a case you would like to try)

To run the GUI you need to have Java installed which you can get free at

https://www.java.com/download/ie_manual.jsp

Also, if you have R (<https://www.r-project.org/>) installed, from the RADDOSE-3D output you will be able to produce 3D representations of the dose distribution in your sample.

Step 3: Find the file RD3D_GUI.EXE and if on a PC click on it. For Linux run it however you usually run executable files. The GUI should open, and you can enter input on 3 tabs: crystal, beam and wedge.



I thank my past and present group and our collaborators,
and acknowledge their huge contribution to the work



Some Garman Group PhD students
and Postdocs 2000-2018:
next generation Crystallographers...



RD Graduate students

James Murray (Imperial College)

Robin Owen (DLS)

Eugenio de la Mora Lugo (IBS, Grenoble)

Oliver Zeldin (Facebook)

Markus Gerstel (DLS)

Helen Ginn (CFEL, Hamburg)

Jonathan Brooks-Bartlett (Spotify)

Charlie Bury (Medicines Discovery Catapult)

RD Postdocs

Karthik Paithankar (U. Frankfurt)

RD Undergraduate Project students

Kathryn Shelley (RABDAM)

Josh Dickerson (RADDPOSE-3D)

RD Collaborators:


Raimond Ravelli, Maastricht.

Martin Weik, IBS

Ian Carmichael, NDRL, USA


John McGeehan, Portsmouth

James Holton, UCSF



The B_{Damage} metric

- There is a strong correlation between mobility and packing density
- Correcting B -factor for packing density enables the distinction of damage from mobility





The B_{Damage} metric

- $B_{\text{Damage}}^{[1]}$ is B -factor corrected for packing density

$$B_{\text{Damage } j} = \frac{B\text{-factor } j}{\frac{1}{n} \sum_{i=1}^n B\text{-factor } i}$$

[1] Gerstel, Deane, Garman (2015) *J Synchrotron Radiat* **22**, 201–212.




Gaussian Beam Profile



Differential irradiation may lead to differential damage:
get data which merge poorly and are a population of substates

13.2 keV, 60(v)×40μm²(h) FWHM, 100(h)×160μm²(v) coll., 5×10^{11} ph/s

Top-Hat profile



Imaged beam
PETRA III,
Bourenkov,
Schneider

13.2 keV, 100(h) x 160 μm^2 (v) coll., 5e11 ph/s