

# CCP4 DLS Workshop 2024

**VERIFICATION:**  
can you prove  
yourself wrong?

**Iracema Caballero**

*icmcri@ibmb.csic.es*

**Arcimboldo team**

Prof. Dr. Isabel Usón

*chango.ibmb.csic.es*


# ARCIMBOLDO: *ab initio* phasing approach

ARCIMBOLDO combines:

- Fragment location with *Phaser* <sup>1</sup>
- Density modification and autotracing with *SHELXE* <sup>2</sup>

Main idea:

- Compose minimal hypotheses from fragments.
- Extend partial structures into full structures with density modification and model building.
- Correct solutions are identified when the minimal starting hypothesis reveals previously unknown structural features, while incorrect ones do not progress into meaningful solutions.



1. A.J. McCoy, R. W. Grosse-Kunstleve, P. D. Adams, et al. "Phaser crystallographic software" (2007) *J.Appl.Cryst*, 40, 658–674  
2. I. Usón and G.M. Sheldrick. "Modes and model building in SHELXE". (2024) *Acta Cryst. D80*, 4–15





# Effect of density modification from a small, accurate fraction of ASU

Phases  
calculated  
from a partial  
solution

Density  
modification

Autotracing of  
new polyalanine  
residues

Iterate: perform  
several cycles



# ARCIMBOLDO programs & modes


## LITE



Single models such as ideal polyalanine  $\alpha$ -helices

- coiled\_coil mode


## SHREDDER



Fragment library from homologs or predicted models

- coiled\_coil mode
- predicted\_model mode

## BORGES



Fragment library representing a fold

## LITE

Sequential search for single models

Signal to noise will be better upon each placement

## SHREDDER & BORGES


Simultaneous search for multiple models of same size

Results are considered jointly



predicted\_model mode in  
ARCIMBOLDO\_SREDDER


# predicted\_model mode in ARCIMBOLDO\_SHREDDER



Medina, A., Jiménez, E., Caballero, I., et al. "Verification: model-free phasing with enhanced predicted models in ARCIMBOLDO\_SHREDDER" (2022). Acta Cryst. D 78(11), 1283–1293


# Model preparation of predicted models

Predicted model




ARCIMBOLDO\_SHREDDER  
performs an **automatic**  
preparation of the input  
template

Pre-processing



- B-factors are set to a common value of 25Å for the main chain and 50Å for the side chains.
- Side chains from predicted models are preserved.
- Hydrogen atoms are removed.


Annotation & Decomposition



- Unstructured areas predicted with low confidence are removed.
- Domains are differentiated through hierarchical decomposition.
- The model is shred into equal-sized spheres based on eLLG.

# Verification

If solution is very partial




**Difficult** phasing is self-validating as wrong solutions impede expansion



Phase combination of partial solutions & expansion

If solution is straightforward




**Straightforward** returning not much beyond the original model



Model-free verification

If coiled coil




**Deceptive** phasing as in coiled-coil structures



Coiled-coil verification


# Model-free verification: eliminating the model bias

If solution is straightforward



In predicted model fragments are now more accurate, but this also introduces model bias

## Model-free verification



Partial solutions are expanded by density modification with SHELXE

One cycle of auto-tracing with SHELXE **omits** the starting model


All consistent traces are combined in reciprocal space with ALIXE, this new map will be iteratively modified and traced

Also used to verify external MR solution

coiled\_coil mode in  
ARCIMBOLDO\_LITE &  
ARCIMBOLDO\_SREDDER

# Coiled coils difficulties for phasing


- Coiled coil structures typically present intrinsic difficulties for phasing:
  - Tendency to aggregate into fibers complicating the formation of well diffracting crystals
  - Large dimensions of their asymmetric units
  - Anisotropic diffraction
  - Highly modulated diffraction patterns
- Phasing of coiled-coil with fragments has been implemented in the AMPLE<sup>1</sup> and CCsolve<sup>2</sup> pipelines.



1. Thomas, J. M. H., Keegan, R. M., Bibby, J., Winn, M. D., Mayans, O. & Rigden, D. J. (2015). IUCrJ, 2, 198-206.

2. Rämisch, S., Lizatović, R. & André, I. (2015). Acta Cryst. D71, 606-614.

# coiled\_coil mode



**ARCIMBOLDO on coiled coils**  
Caballero, I., Sammito, M., Millán, C., et al.  
Acta Cryst. **2018**; D74, 194–204

Up to 3Å

**ARCIMBOLDO at low resolution: Verification for coiled coils and globular proteins**  
Caballero I, Castellví A, Triviño J, et al.  
Protein Sci **2024**; 33(9):e5136.

Up to 4Å




| PROBLEMS  | SOLUTIONS  |
|---|--|
| Difficulties to differentiate genuine intermolecular tNCS from Patterson artefacts. | Disable the placement of pairs of tNCS related helices.                            |
| Overlapping solutions with high FOMs.   | Activate the <i>Phaser</i> packing filter during the translation search.           |
| The geometry of helical trace can degrade.  | SHELXE with helical sliding that anchors ahead of the trace, improving the phases. |
| Placement of helices in the correct position but in reversed direction.             | Generate and test reversed helices in the same positions (only in LITE).           |
| Wrong solutions can have high figures of merit.                                     | Final verification step.   |





# Resolution limits

Resolution in ARCIMBOLDO is **auto-limited** in order to prevent false positives




# Verification step - prove ourselves wrong

To verify the most promising solution, its original substructure will be perturbed by a random translation and by reversing the direction of the helices




In LITE & SHREDDER




← Only in LITE

# Verification in ARCIMBOLDO\_LITE

The extensions are compared for evidence of discrimination between the randomly translated solution or groups of consistent solutions




The best solution is **discriminated** from the random solution or the perturbations develop into a group of structurally **equivalent solutions**




The best solution is **not discriminated** from the random solution or the perturbations develops into structurally **different solutions**

# Verification in ARCIMBOLDO\_SHREDDER

The extensions are compared for evidence of discrimination between the best solution and the randomly translated solution





The best solution is discriminated from the random solution




The best solution is not discriminated from the random solution

# Verification as a key tool to eliminate false positives

Verification ruled out 7 false positives in the case of ARCIMBOLDO\_LITE and 4 false positives in the case of ARCIMBOLDO\_SHREDDER




# Choices for phasing coiled coils



ARCIMBOLDO\_LITE in coiled\_coil mode

## Choice for searching single helices

It performs a sequential search for helices of a given length with specific features to avoid coiled-coil pitfalls.



ARCIMBOLDO\_SHREDDER

in coiled\_coil mode & predicted\_model mode

## Choice for searching multimeric predictions

It extracts compact fragments from the starting model and gives them degrees of freedom to improve the models.

Low resolution verification for  
globular proteins



# Structure approximation by a Fourier series

- Structure determination by crystallographic diffraction can be expressed through the Fourier series.
- This animation uses a line drawing of a hummingbird to define a periodic path, where the coordinates in the complex plane are a function of time. The path can be approximated by a Fourier series equation:

$$F(t) = \sum coef_n \cdot \exp(i \cdot 2\pi \cdot n \cdot t)$$

- The Fourier series approximation is only exact with infinite terms.
- With fewer terms (low resolution), we not only lose detail but may also introduce artifacts that do not exist in the original structure.

<http://chango.ibmb.csic.es/colibri>







# LLG-based fragment verification from density modified maps

1. Constrained map interpretation: Use SHELXE density-modified maps from partial, reliable solutions to build hypotheses consistent with prior knowledge.
1. Comparative scoring of alternative hypotheses: Evaluate fragments from remote homologs or secondary structure predictions:
  - Fragments were rigid-body refined, allowing internal degrees of freedom using *Gimble* in *Phaser*.
  - Fragments were scored by their  $LLG_{\text{contribution}}$  to the overall value, estimated from LLG change upon fragment omission.
  - Relative differences in the  $LLG_{\text{contribution}}$  of each fragment will guide the choice between different possibilities: size, orientation, secondary structure or side-chain inclusion.


# Multidomain protein VSR1 at 3.5Å



The complex state of the PA domain 4txj rendered a partial solution of VSR1 by MR with *Phaser*.

No models were available for the Central domain, but has two remote homologs with 15% identity.


Both structures share a roughly parallel four-helix bundle orientation.




Staphylococcus aureus DsbA (PDBid 3bci)  
Aspergillus fumigatus (PDBid 4k90)

# Multidomain protein VSR1 at 3.5Å

Density-modified  
map calculated  
using the  
coordinates of  
the PA domain




The four-helix bundle was placed in the Central  
domain region and scored on LLG<sub>contribution</sub>




Optimization of the blue helix orientation  
increases its LLG<sub>contribution</sub>

# Multidomain protein VSR1 at 3.5Å




Hypothesis of four helices of the VSR1 superimposed to DsbA (15% identity), the green helix in the central domain corresponds to a strand of DsbA.


We included the central domain fragments in the map calculation



The green strand was expanded to a sheet of 3 strands, and adding side chains improved their LLG<sub>contribution</sub>.



The green helix was replaced by a strand and its LLG<sub>contribution</sub> increases.




Density-modified map at the EGF region generated solely with the PA domain vs. using also fragments from the Central domain. **Emergence of new features in unmodeled map regions confirmed partial model correctness.**




Tutorial

# ARCIMBOLDO Interfaces



## CCP4i




## CCP4i2




## CCP4 online (Ville Uski)



## CCP4 Cloud



## XDSGUI (Kay Diederichs)



# Setting the job in CCP4 online

## Programs

Note: You must have a CCP4 licence to run these programs.

|            |  |
|------------|--|
| Balbes     | An automated Molecular Replacement (MR) pipeline - Balbes integrates into one system all the components necessary for solving a crystal structure by Molecular Replacement   |
| MrBUMP     | An automated Molecular Replacement (MR) pipeline - Given a target sequence and experimental structure factors, it will search for homologous structures, create a set of suitable search models from the template structures, do molecular replacement, and test the solutions with some rounds of restrained refinement.<br>New: MrBUMP now searches the EBI-AFDB AlphaFold database for potential search models in addition to the PDB. This increases the pool of potential search models by 350000 with more to come in the near future. |
| Zanuda     | Space group and crystallographic origin validation   |
| AMPLE      | Automated ab initio search model generation for molecular replacement.   |
| SHELX      | Automated SHELXC/D/E structure solution pipeline for fast routine experimental phasing. Accepts data in XDS, Scalepack, SHELX hkl or mtz formats and outputs phases and a poly-Ala trace. If a protein sequence is provided, BUCCANEER and REFMAC complete the structure.  |
| CRANK2     | Automated structure solution pipeline for experimental phasing using maximum likelihood methods.   |
| MoRDa      | MoRDa is a pipeline for molecular replacement protein structure solution based on its own domain database. Models relevant to the target sequence are further adjusted before molecular replacement search.  |
| SIMBAD     | Sequence-independent molecular replacement, good for identifying if your crystal contains a contaminant protein. SIMBAD can also do full search of homologous structures in difficult-to-solve novel target cases, but this functionality is not yet available through CCP4-Online.  |
| ARCIMBOLDO | Crystallographic Ab Initio protein solution below atomic resolution  |

### ARCIMBOLDO

[CLICK HERE TO START A NEW ARCIMBOLDO RUN](#)

The table below shows your current ARCIMBOLDO runs - the results will be stored on the server for at least 30 days

| id | Job Title | Delete/Stop | Date | Status | View Results |
|----|-----------|-------------|------|--------|--------------|
|----|-----------|-------------|------|--------|--------------|

# Setting the job in CCP4 online

## New ARCIMBOLDO Run

ARCIMBOLDO combines fragment search with PHASER and density modification and autotracing with SHELXE, at resolutions of 2.5 Å or better (4 Å for coiled-coil mode).

ARCIMBOLDO\_LITE uses single models, typically polyalanine helices;

ARCIMBOLDO\_BORGES uses precomputed libraries of local folds;

ARCIMBOLDO\_SHREDDER uses fragments from a distant homolog or predicted models.

### Data description

|                              |                         |  |
|------------------------------|-------------------------|--|
| run ARCIMBOLDO:              | SHREDDER                | ARCIMBOLDO SHREDDER selection                    |
| Job title (optional):        | 8b8d                    |  |
| Input data (.mtz):           | Browse... 8b8d.mtz      | Upload mtz                                       |
| Input type:                  | amplitudes              |  |
| Column labels:               | SIGF FP F SIGFP         | Select labels                                    |
| Target Sequence:             |                         |  |
| Asymmetric unit description: | Molecular weight: 51680 | (please give either this or the target sequence) |
|                              | Number of components: 1 | Molecular_weight = 51680 Da                      |

number\_of\_component = 1

## ARCIMBOLDO SHREDDER

|  |                            |                               |
|--|----------------------------|-------------------------------|
| Model:                                 | Browse... 8b8d_model_0.pdb | r.m.s.d.: 0.8 Å               |
| Convert to polyalanine (optional):     |                            |                               |
| Make all B-factors equal (optional):   |                            |                               |
| Alphafold or RoseTTaFold model:        | True                       | Activate predicted_model mode |
| Shredder mode:                         | spherical                  |                               |
| Fragment size (optional):              |                            |                               |
| Maintain coil in the model (optional): |                            |                               |
| Perform gyre refinement (optional):    |                            |                               |
| Perform gimble refinement (optional):  |                            |                               |
| Perform LLG-guided pruning (optional): |                            |                               |
| Combine phases with alixie (optional): |                            |                               |

Upload AF model

Activate predicted\_model mode

### Advanced options

|                   |             |            |
|-------------------|-------------|------------|
| Add keyword-value |             |            |
| Keyword           | coiled_coil | Value True |
|                   |             | remove     |

Activate coiled\_coil mode

|        |          |         |
|--------|----------|---------|
| Submit | Validate | Restart |
|--------|----------|---------|

Click Validate

Then click submit



# Running the job: HTML results

The table below shows your current ARCIMBOLDO runs - the results will be stored on the server for at least 30 days

| id         | Job Title | Delete/Stop            | Date                | Status        | View Results            |
|------------|-----------|------------------------|---------------------|---------------|-------------------------|
| 1552843156 | 8b8d      | <a href="#">delete</a> | 12 Nov 2024 [11:14] | process ended | <a href="#">results</a> |

## ARCIMBOLDO

PROCESS 1552843156 HAS ENDED

### ARCIMBOLDO-SHREDDER

Summary

Instruction file

Log file

Citations

Support

#### Summary

NAME: 8b8d  
SPACEGROUP: P 21 21 2  
CELL DIMENSIONS: 81.82, 82.50, 89.06, 90.00, 90.00, 90.00  
RESOLUTION: 2.40  
NUMBER OF UNIQUE REFLECTIONS: 19316.00

[Click here](#) to follow the progress with the ARCIMBOLDO\_BORGES

It seems you have a good solution!  
Here you can find the best [solution](#) and [map](#) for further refinement.

#### Instruction file

Input file at [/ccp4online/ARCIMBOLDO/iracaballero/1552843156/input/instructions.bor](#). Click [here](#) to download it

```
[CONNECTION]
distribute_computing = multiprocessing
remote_frontend_passkey = ~/.ssh/id_rsa
setup_bor_path =
```

# Running the job: HTML results

## ARCIMBOLDO-BORGES

### Search and expansion

| #Cluster | Rotation and Model Refinement |         |          |            |             | Translation Function |         |          |            |             | Packing |         |          |            |             |       |
|----------|-------------------------------|---------|----------|------------|-------------|----------------------|---------|----------|------------|-------------|---------|---------|----------|------------|-------------|-------|
|          | #Rots.                        | Top LLG | Mean LLG | Top Zscore | Mean Zscore | #Trans.              | Top LLG | Mean LLG | Top Zscore | Mean Zscore | #Sol.   | Top LLG | Mean LLG | Top Zscore | Mean Zscore | #Sol. |
| 1        | 459                           | 243.00  | 133.97   | 4.16       | 2.83        | 100                  | 347.04  | 188.48   | 13.97      | 6.40        | 100     | 347.04  | 188.48   | 13.97      | 6.40        | 100   |
| 5        | 446                           | 229.00  | 131.85   | 4.23       | 2.77        | 100                  | 252.01  | 156.42   | 8.64       | 4.78        | 100     | 252.01  | 156.42   | 8.64       | 4.78        | 100   |
| 0        | 447                           | 244.00  | 134.25   | 4.18       | 2.79        | 100                  | 336.31  | 181.62   | 12.68      | 6.08        | 100     | 336.31  | 181.62   | 12.68      | 6.08        | 100   |
| 6        | 438                           | 242.00  | 132.01   | 4.23       | 2.75        | 100                  | 349.79  | 156.82   | 13.40      | 4.75        | 100     | 349.79  | 156.82   | 13.40      | 4.75        | 100   |

Show AllHide Not Relevant

### Backtracking

The current best solution is: frag41\_0\_0.pdb with FINALCC: 54.48 and n. residues traced 367  
file is: /ccp4online/ARCIMBOLDO/iracaballero/1552843156/ARCIMBOLDO\_BORGES/.11\_EXP\_alixe/1/5/0/frag41\_0\_0\_rbr\_1\_ref.pdb

#### • REFINEMENT ROTATION AND MODEL

#### • EXPANSION


Cycle 5  
Final CC: **54.48%** N. Residues Traced: **367.00**

### Verification


Verification has determined that the structure is solved.  
The best solution is clearly distinguished from a random one.

Solved!

Verification step



PDB entry: 8B8D  
Resolution: 2.4Å  
4 helices of 115 residues



# SUMMARY - ARCIMBOLDO

- ARCIMBOLDO\_SHREDDER has a **predicted\_model mode**:
  - No need to preprocess models.
  - If solution is straightforward, it will systematically eliminate model bias.
  - This feature can be run as standalone from a MR solution.
- The **coiled\_coil mode** is available in ARCIMBOLDO\_LITE and ARCIMBOLDO\_SHREDDER:
  - Works at resolution up to 4Å.
  - Use ARCIMBOLDO\_LITE with ideal helices and ARCIMBOLDO\_SHREDDER with multimeric predictions.
  - It includes a verification step, key when going to lower resolutions.

# SUMMARY - VERIFICATION

The core idea of ARCIMBOLDO is that correct solutions are identified when the minimal starting hypothesis reveals previously unknown structural features, while incorrect ones do not progress into meaningful solutions.

- **Predicted models:** the starting hypothesis leaves little to infer as fragments are now more accurate but this also introduce model bias → model-free verification that will systematically eliminate model bias.
- **Coiled coils:** figures of merit become unreliable → dedicated verification step to disprove the final solution against competing alternatives.
- **Low resolution in globular proteins:** the emergence of new features in unmodeled map regions confirmed partial model correctness.

# ACKNOWLEDGEMENTS



ARCIMBOLDO TEAM: **Isabel Usón**, Josep Triviño, Elisabet Jiménez, Mari Carmen Martín & Iain Visa

~~Old~~ Always members of the ARCIMBOLDO TEAM: Claudia Millán, Massimo Sammito, Albert Castellvi, Ana Medina, Nicolas Soler & Rafael Borges

PHASER group: Airlie McCoy & Randy Read

Kay Diederichs & George Sheldrick

